




Pride® Lift Chairs
Experience the Luxury

Protect your
investment with
the Pride® Lift
Chair Care Kit



The Specialty Collection from Pride® offers a wide selection of high-quality lift chairs. Many chairs in the collection have distinctive options and special features ranging from enhanced dimensions to infinite positioning and much more. The Specialty Collection also contains special value chairs that deliver unique options at an affordable price.

SPECIALTY COLLECTION

www.pridemobility.com

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

Pride®
Mobility Products Corp.

SPECIALTY COLLECTION MODELS

SPECIAL FEATURES

							
Small	Large	Medium	Extra Tall	Large	Medium	Medium	Medium
LL-770S	LL-770L	LL-570T⁹	LL-575⁹	LL-571^{1,3,9}	LL-805⁴	GL-310	LC-900^{6,7,8}
Infinite-Position Trendelenburg Chaise Lounger 5'4" - Below	Infinite-Position Trendelenburg Chaise Lounger 5'5" - 6'0"	3-Position Chaise Lounger w/ Transfer Arms 5'5" - 6'0"	3-Position Chaise Lounger w/ Footrest Extension 6'1" - 6'6"	3-Position Chaise Lounger w/ Flip-open Arms 5'5" - 6'0"	2-Position Wall Hugger Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'2" - 5'10"	3-Position Chaise Lounger 5'5" - 6'0"

**SPECIALTY
COLLECTION**

SPECIAL VALUE





				
Super Petite	Medium	Small	Medium	Large
LC-102	LC-105	LC-125S	LC-125M	LC-125L
3-Position Chaise Lounger 4'8" - 5'4"	3-Position Chaise Lounger 5'4" - 6'0"	Infinite-Position Trendelenburg Chaise Lounger 5'4" - Below	Infinite-Position Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite-Position Trendelenburg Chaise Lounger 5'5" - 6'0"

SPECIAL FEATURES FABRIC SELECTION

STANDARD FABRICS				
				
MARINE	WINE	MOSS	WHEAT	TAUPE
OPTIONAL FABRICS				
				
RASPBERRY	HERB	ALMOND	TRUFFLE	
				
BRASS	ROSE TWIST	OATMEAL	SLATE	
LEXIS VINYL				
				
BLACK	MUSHROOM	BURGUNDY	NEW CHESTNUT	ROYAL BLUE
ULTRAFABRICS® ULTRALEATHER™				
				
BUFF	FUDGE	GARNET	NILE	PECAN
CRYPTON® SUPER FABRICS				
				
OYSTER	PICKLE	CLOVER	MALT	GINGER
SMART FABRICS**				
				
VANILLA BEAN	RED PEPPER	CELERY	CREAM	
NEW OUTLAST® TECHNOLOGY				
				
CHAMPAGNE	NUTMEG			

SPECIAL VALUE FABRIC SELECTION



LC-105 FABRIC SELECTION

			
COCOA	FERN	SANDAL	SKY

LC-125 SERIES FABRIC SELECTION

				
NEW CHESTNUT VINYL	BLACK VINYL	SAND	OCEAN	TEABERRY

LC-102 FABRIC SELECTION

	
COCOA	SANDAL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

**5-year warranty.

SPECIFICATIONS

SPECIAL FEATURES

	LL-770S	LL-770L	LL-570T	LL-575	LL-571	LL-805	GL-310
Seat-To-Floor:	19"	21"	20"	21"	21"	19"	20"
Seat Depth:	19.5"	21"	20"	22"	22"	20"	20"
Seat Width:	19"	22"	20"	22"	22"	20"	21"
Top of Back to Seat:	24"	28"	29.5"	31"	27"	27"	24.5"
Distance From Wall:	28"	30"	20.5"	24"	16"	4"	18"
Overall Chair Width:	34"	37"	35.5"	41.5"	38"	35"	30"
Chair Weight	129 lbs.	136 lbs.	128 lbs.	138 lbs.	133 lbs.	118 lbs.	115 lbs.

SPECIAL FEATURES

SPECIAL VALUE

	LC-900	LC-102	LC-105	LC-125S	LC-125M	LC-125L
Seat-To-Floor:	19"	16"	18"	18.5"	19"	19"
Seat Depth:	20"	17"	19"	18"	20"	19"
Seat Width:	22"	16"	20"	20"	20"	22"
Top of Back to Seat:	22"	23.5"	27.5"	23"	28.5"	28.5"
Distance From Wall:	17"	16"	18"	19"	24"	24"
Overall Chair Width:	34"	25.5"	32.75"	31"	31.5"	35.5"
Chair Weight	120 lbs. ⁵	96 lbs.	96 lbs.	126 lbs.	136 lbs.	140 lbs.

BENEFITS

- Engineered furniture grade laminate/hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head & arm covers (head cover only on GL-310, LL-770 Series and LL-900; no head or arm covers on LC-105 or LC-125 Series)
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets (not available on LC-900)
- Removable back - no tools required (not available on LC-125 Series)
- 375 lbs. weight capacity (LL-805 supports 300 lbs.; LC-102 supports 250 lbs.; LC105 supports 325 lbs.)
- Integrated battery backup in case of power failure²
- Transformer features self-diagnostic electronics with 12' cord
- 4-year limited warranty on LL-805
- 7-year prorated warranty (2-year limited warranty on LC-102 and LC-125 Series; 1-year limited warranty on LC-105)

OPTIONS*

Additional lead time may apply for some options.

- Heat & Massage
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Bed option w/mattress, fitted sheet, & storage bag (LL-770L only)
- Non-skid leg levelers
- Leather/vinyl wipes
- Media wedge (LC-900 only)

1. Standard with compact table that stores conveniently in chair arm.
2. Batteries not included.
3. Removable, fully padded chaise.
4. Non-removable (adjustable) chaise pad.
5. Media wedge weight: 38 lbs.
6. Only available in Lexis Vinyl.
7. Media wedge dimensions: 12"W x 26.5"H x 32"L.
8. No side pockets.
9. Fiber fill can be added or removed for custom comfort adjustment.

*Additional charges apply for all options. Not all options are available on every model. Please consult order form for compatibility.

Look for the Quick Ship Indicator for the shortest lead times available.