


Arlo Rotating Bed

Key Features


- Multifunctional Design
- Silent and Durable Motor
- One-Touch 90° Rotation Bed
- Head (90°) and Foot (35°) Elevation
- Air Fiber Breathable Mattress
- Sideboards and Safety Belt For Protection
- US FDA Registered For Bed and Mattress
- 450lbs Weight Capacity

2 Bed Models for Selection

Left Rotation | Right Rotation



All-In-One Features



Product	Volume	Product Dimension	Package Dimension	N.W/G.W	
Arlo	47.3 ft ³	91x39x40in	84.6x37.4x16.3in	211lbs	218lbs
			67.3x38.9x13in	101lbs	108lbs