





Ventura 4 DLX Midsize Scooter




INSET A



INSET B



INSET C



INSET D



- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- 33AH batteries for extended range
- Ergonomic throttle control
- Height-adjustable swivel seat adjusted forward and backward features fold-down backrest and adjustable headrest
- Front and rear suspension with Anti-tip wheels
- Programmable controller
- Available with 18" or 20" width Captain's Seat
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, rearview mirrors and turn signals/hazard lights
- Easy free wheel operation

PRODUCTS	Description
VENTURA418CS	1/cs
VENTURA420CS	1/cs

SPECIFICATIONS	Description
HCPCS	N/A
Weight Capacity	400 lbs.
Top Speed	6 mph
Maximum Range	16 miles
Turning Radius	54"
Climbing Angle	7.5°
Ground Clearance	2.5"
Dimensions	48" (L) x 25" (W)
Floor-to-Seat Height	22" - 24"
Seat Type	Captain's
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Motor	24V x 350W x 5100 rpm
Batteries (pair)	12V x 33AH
Charger	4A Offboard
Brakes	Electromagnetic
Front Wheels	10" x 3" Flat-free
Rear Wheels	10" x 3" Flat-free
Base Weight	99 lbs.

SPECIFICATIONS	Description
Battery Weight (pair)	46 lbs.
Seat Weight	40 lbs.
Heaviest Piece Weight	70 lbs.
Total Weight	185 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Batteries	12 Months

ACCESSORIES	Description
ACCKIT-M	1/cs
SF8010	1/ea
SF8020	1/ea
SF8030	1/ea
SF8050	1/ea

PARTS	Description
VENT-01	1/ea
VENT-02	1/ea
VENT-03	1/ea
VENT-04	1/ea
VENT-05	1/ea
VENT-06	1/ea
VENT-07	1/ea
VENT-08	1/ea
VENT-09	1/pr
VENT-10	1/ea
VENT-11	1/ea
VENT-12	1/ea
VENT-13	1/ea
VENT-14	1/ea
VENT-15	1/ea
VENT-44	1/ea
VENT-17	1/ea
VENT-16	1/ea
VENT-18	1/ea
VENT-19	1/ea
VENT-20	1/ea
VENT-21	1/ea
VENT-22	1/ea
VENT-36	1/ea
VENT-23	1/ea
VENT-26	1/ea
VENT-27	1/ea
VENT-24	1/ea
VENT-25	1/pr
VENT-28	1/ea
VENT-29	1/ea
VENT-30	1/ea
VENT-31	1/ea
VENT-32	1/ea
VENT-33	1/ea
VENT-34	1/ea
VENT-35	1/ea
VENT-41	
VENT-42	