

INVACARE®

# Veranda™

Wheelchair



Yes, you can.®

The **Invacare® Veranda™ Wheelchair** provides comfort and mobility in an economical package. It comes complete with front riggings packaged in the carton with the chair. The Veranda Wheelchair combines style and durability at an affordable price.



## Features

- Padded nylon upholstery
- Comfortable padded arms
- Stylish powder-coated frame
- Includes front riggings
  - Swingaway Footrests match chair color
  - Elevating Legrests are chrome
- Carry pocket on back
- Large step-tube end cap
- Dual axle positions allow variation of seat-to-floor heights

**Invacare Corporation**  
www.invacare.com

## USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

## Canada

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

© 2010

Invacare Corporation.

All rights reserved.

Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.

Form No. 07-064  
Rev. 10/10  
100947



Pocket on back.



Step-tube.



Large step-tube end cap.

## Specifications

|   |                            |
|---|----------------------------|
| <b>HCPCS code</b>                       | K0001/K0002                |
| <b>Seat sizes</b>                       |                            |
| Width                                   | 16", 18", 20"              |
| Depth                                   | 16"                        |
| <b>Seat-to-floor height</b>             |                            |
| Removable arm                           | 17.5", 19.5"               |
| Permanent arm                           | 17.5", 19.5"               |
| <b>Back height</b>                      | 18"                        |
| <b>Front riggings</b>                   |                            |
| Hemi swingaway footrest                 | 16.5" - 19.5" knee-to-heel |
| Hemi swingaway padded elevating legrest | 17.5" - 21" knee-to-heel   |
| <b>Arm styles</b>                       |                            |
| Permanent full-length, fixed height arm |                            |
| Removable desk-length, fixed height arm |                            |
| <b>Weight capacity</b>                  | 250 lb.                    |
| <b>Product weight</b>                   | 37 lb.                     |
| <b>Limited warranty</b>                 | 1 year                     |

