


PRELIMINARY *viva*LIFT!™



Shown in Merino
Fawn Fabric



Shown in Badlands
Walnut Fabric

*viva*LIFT!™ POWER RECLINERS

By Pride® Mobility Products Corp.

Breakthrough Technology!

Power Lumber for personalized ergonomic comfort

Full-Width Power Headrest to lift head,
neck and shoulders

Lithium Battery Back-Up is an industry first!
Abundant power lasts for days and replaces
ordinary one-time-use backups

Fashion Forward Styles set these chairs
apart from all others!

Elevate Your Lifestyle!

(US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com

Pride®
Mobility Products Corp.
— LIVE YOUR BEST! —

VivaLIFT!™

POWER RECLINERS

By Pride® Mobility Products Corp.

Elegance



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

Urbana



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

Perfecta



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

Metro



STANDARD FEATURES:

- POWER HEADREST
- POWER LUMBAR
- USB REMOTE
- LITHIUM BATTERY BACKUP
- ALL CLOSE AND LIFT
- SPRING LOADED OTTOMAN
- INFINITE LAY FLAT

Model	PLR975M	PLR965M	PLR945M	PLR925M
Position	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'	Relaxed Reclined Infinite Lay Flat Upward Lift 5'4" - 6'
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	20"	19"	19"
Seat Depth	20"	20"	20"	21"
Seat Width	22"	21"	22"	21"
Top of Back to Seat	27"	25"	25"	26 1/2"
Back Style	Blown fiber, Pillow back	Blown fiber, Split back	Blown fiber, Split button back	Blown fiber, Split back
Fully Padded Chaise	Yes	Yes	Yes	Yes
Footrest Extension	Standard	Standard	Standard	Standard
Distance From Wall	14"	14"	15"	14"
Height	43"	42"	43"	41"
Overall Width	35 1/2"	33"	31"	33"
Chair Weight	152 lbs.	145 lbs.	142 lbs.	142 lbs.
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

STANDARD FABRIC SELECTION GUIDE



ICON KEY
 = Indicates shortest available lead time

LC-VIVA-SS-09-28-17

(US) 800-800-8586 | (CAN) 888-570-1113
 pridemobility.com | prideliftchairs.com

Pride
 Mobility Products Corp.
 — LIVE YOUR BEST —