

## Vivo Features

VIVO

- **10° Fixed Tilt**  
*Promotes an upright position for individuals with good trunk strength and upper body stability*
- **95° – 120° of Adjustable Recline**  
*Opens the seat-to-back angle for better digestion, improved respiration, and pressure relief*
- **Adjustable Strap Back**  
*Provides contoured positioning*
- **Zippered Pockets in Upholstery**  
*Allow additional foam padding to be added for firmer support*
- **Adjustable Strap Seat Supports**  
*Adjust position of hips to provide pressure relief*
- **Ultra-Lightweight Wheelchair**  
*Is easy to lift and carry*
- **Compact Folding Frame**  
*Easily stows in a car, bus or airplane – fold and go!*
- **Comprehensive Positioning Options**  
*Provide comfort and proper positioning*
- **Multiple Frame Adjustments**  
*Extend the life of the chair for a growing child*



Positioning options and accessories sold separately.

# Vivo Specifications

## Vivo Standard Features

- 2-Point Positioning Belt
- Height Adjustable Swing-Away Footplates
- Closure Strap
- Breathable Seat Back Fabric
- Toggle Wheel Locks
- 6" (15cm) Front – 8" (20cm) Rear Tires
- Pac Cloth Fabric

## Options and Accessories

### Head Positioning

Headrest Extension  
Padded Neck Support

### Pelvis Positioning

Height Adjustable Armrest  
Foam Wedge

### Accessories

Padded Tray  
Canopy  
Grab Bar  
Wheelchair Travel Bag

### Trunk Positioning


Anterior Trunk Support Vest  
H-Harness with Padded Covers

### Foot and Leg Positioning

Foot Positioners  
Swing-Away Hip Abductor  
Calf Panel

## Upholstery Colors

Removable,  
Washable Pac Cloth Fabric



## Vivo Specifications

### Vivo 12

Seat Width: 12" / 30 cm  
Seat Depth: 9.8" - 11.8" / 25 - 30 cm  
Chair Length: 31.5" / 80 cm  
Chair Width: 19.7" / 50 cm  
Chair Height: 36.6" / 93 cm  
Seat Back Height: 21.65" / 55 cm  
Seat-to-Floor Angle: 10°  
Seat Back Angle: 95° - 120°  
Footplate Height: 7.9" - 11" / 20 cm - 28 cm  
Weight Capacity: 44 lbs / 20 kg  
Weight of Chair: 14.3 lbs / 6.5 kg

### Folded Dimension

Folded Length: 38.2" / 97 cm  
Folded Width: 14.6" / 37 cm  
Folded Height: 17" / 43 cm

### Vivo 14

Seat Width: 14" / 35.6 cm  
Seat Depth: 11.8" - 14" / 27.5 - 35 cm  
Chair Length: 39" / 99 cm  
Chair Width: 21.65" / 55 cm  
Chair Height: 41.2" / 104.6 cm  
Seat Back Height: 26" / 66 cm  
Seat-to-Floor Angle: 10°  
Seat Back Angle: 95° - 120°  
Footplate Height: 9.8" - 14" / 25 - 35.6 cm  
Weight Capacity: 88.2 lbs / 40 kg  
Weight of Chair: 17.2" / 7.8 kg

### Folded Dimension

Folded Length: 46.6 lbs / 118.4 cm  
Folded Width: 11.4" / 29 cm  
Folded Height: 17.4" / 44.3 cm



U.S. Patent 4,118,065  
This product may be covered by the  
following U.S. patent - 6,099,075  
Additional patents pending.

visit us at [www.convaid.com](http://www.convaid.com)