



SURRENDER NOTHING.

LIVE LIFE AT
FULL THROTTLE,



XENON²

SFC (Retired) Mark H.
U.S. Army 3rd Special Forces Group

XENON²

LIGHTWEIGHT,
HIGH-PERFORMANCE,
FOLDING SERIES

WHETHER YOU'RE LOOKING FOR
THE VERY BEST IN DRIVING DYNAMICS
OR THE MOST COMPACT FOLDING
DIMENSIONS, THERE'S A XENON² JUST FOR YOU.




LIGHTEST FOLDING WHEELCHAIR ...IN ITS CLASS*

How do you achieve a folding wheelchair that weighs as little as 20 pounds? By using premium frame materials like 7000 series aerospace aluminum and an innovative process that's exclusive only to QUICKIE® – ShapeLoc Technology.

UNRIVALED FOLDING PERFORMANCE

At the core of the XENON² is a unique cross-brace that fits so neatly under the seat it's barely noticeable – giving a very minimalistic, open-frame feel that you'd usually associate with a rigid wheelchair.



* Class is defined as: aluminum folding wheelchairs (North America).

AVAILABLE IN 3 MODELS:



XENON² FIXED FRONT
THE LIGHTEST



XENON² HYBRID (DUAL TUBE)
THE STRONGEST



XENON² SWING-AWAY
THE ALLROUNDER



Sunrise Medical (USA) LLC, Fresno, CA 93727
MK-129596 Rev.A ©10.2017

800.333.4000 www.SunriseMedical.com

