



Invacare® Xtra™ Wheelchair



Yes, you can.™



Yes, you can.™

The Invacare® Xtra™ chair is designed to offer maximum performance, increased stability, and everything in between. This range of handling characteristics is achieved through the combination of the extended rear frame and an axle plate that delivers a remarkable 66 different wheel positions. The Xtra also features a superlite frame which is split to facilitate growth. All this at a price that brings the high end of the folding wheelchair market within reach.

Highlights

- Box rear frame
- Two-piece side frames
- More features than any other folding chair in its price category
- Choice of two-point or swingaway arms
- Choice of seven different front riggings

Specifications

Dimensions		
Overall		
Width open	8.25" plus seat width*	
Front Seat-to-Floor		
Fork type	6" caster	8" caster
6" upper mounting holes	17", 17.75", 18.5"	n/a
6" lower mounting holes	18.25", 19", 19.75"	n/a
8" 2-holes	n/a	18.75", 20", 21.25"
Rear Seat-to-Floor		
Wheel Size (range adjusts in 7/8" increments)		
20"	14.75"-17.75"	
22"	15.75"-18.75"	
24"	16.5"-19.75"	
Knee-to-Heel** (range adjusts in 1" increments)		
70° tapered	15.75"-19.75"	
70°	13.5"-18.5"	
70° MFX	8.5"-14.75"	
Elevating	15.25"-21.25"	
90°	5.25"-10.25"	
60°	14.75"-18.5"	
60° MFX	11"-15"	
Arm Heights*** (range adjusts in 1" increments)		
Adjustable-height dual-point	10"-14" spacer saver only	
Fixed-height dual-point	10" spacer saver only	
Swingaway half-arm	8"-10" add 1" with upholstered arm pads	
Cane Options		
Adjustable-height	Sportster, 10° with push handles, Straight with push handles	
Fixed-height	10° with push handles 20" height	

*Measurement taken from outside of handrim to outside of handrim with 0° camber.

**Measurement taken from top of seat rail to back of footrest.

***Measurement taken from top of seat rail to top of armrest.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

5970 Chedworth Way
Mississauga, Ontario
L5R 3T9 Canada
(800) 668-5324

Invacare is a registered trademark and Yes, you can. and Xtra are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2002
Invacare Corporation

Form No. 00-50
Rev. 05/02



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."