

ZooMe Flex Folding Travel Scooter



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Ultra-lightweight for superior portability
- Quick and Easy "Flex" folding design for transport and storage in one easy motion
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests and fold down back
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position
- Anti-tip wheels for additional safety

PRODUCTS

Description

FLEX	ZooMe Flex Folding Travel Scooter, Blue
------	---

SPECIFICATIONS

Description

Anti-Tip Wheels	2" x .5"
Base Weight	53.6 lbs.
Batteries (pair)	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Battery Weight (pair)	5 lbs.
Brakes	Electromagnetic
Climbing Angle	6°
Controller	PG S Drive 45A
Estimated Range	13 miles
Freewheel Mode	No
Front Wheels	6" x 1.5"
Ground Clearance	2"
HC-FEET	4" (L) x 3" (W) x 0.5" (H)
HC-HAND	5" (L) x 3" (W) x 1.25" (H)
Heaviest Piece Weight	53.6 lbs
Motor	24V x 120W
Overall Dimensions	37" (L) x 18" (W) x 36" (H)
Rear Wheels	7" x 2.35"
Seat Dimensions	16.75" (W) x 15" (D)
Top Speed	4 mph
Total Weight (excl. hand control)	58.6 lbs.
Turning Radius	53"

SPECIFICATIONS**Description**

Warranty	24 months on Electronics, 12 months on Batteries, Limited Lifetime on Frame
Weight Capacity	275 lbs
Weight Capacity (13046 13046-HD)	275 lbs.

PARTS**Description**

FLEX-01	Upper Panel Assembly,Blue,ZooMeFlex,1/ea, EA
FLEX-02	Throttle Assembly, ZooMe Flex, 1/ea, EA
FLEX-03	Gas Spring Assembly, ZooMe Flex, 1/ea, EA
FLEX-04	FrontReflectorAssmbly,Blue,ZooMeFlex,1ea, EA
FLEX-05	Throttle Potentiometer, ZooMe Flex, 1/e, EA
FLEX-06	Quick Release Clamp, ZooMe Flex, 1/ea, EA
FLEX-07	Front Wheel Assembly, ZooMe Flex, 1/ea, EA
FLEX-08	Steering Knuckle, Left, ZooMe Flex, 1/e, EA
FLEX-09	Steering Knuckle,Right, ZooMe Flex, 1/e, EA
FLEX-10	Hub Cap, ZooMe Flex, 1/ea, EA
FLEX-11	Wheel / Stem Bearings, ZooMe Flex, 1/pr, EA
FLEX-12	Front Cover, Blue, ZooMe Flex, 1/ea, EA
FLEX-13	Tie Rod W/Ends, ZooMe Flex, 1/ea, EA
FLEX-14	Handle Connector, ZooMe Flex, 1/ea, EA
FLEX-15	Steering Connector, ZooMe Flex, 1/ea, EA
FLEX-16	BatteryCase W/O Batteries,ZooMe Flex,1e, EA
FLEX-17	Front Pedal Cover, Blue,ZooMe Flex, 1/ea, EA
FLEX-18	RotatingGearAssembly,Front,ZooMeFlex,1e, EA
FLEX-20	ReleaseCableW/Hrdwr&Lever,ZooMeFlex,1/ea, EA
FLEX-19	Support Rod, Front, ZooMe Flex, 1/ea, EA