

# aXs Lift



The aXs lift is a perfect accessory for any backyard pool, where its low profile and attractive appearance make swimming pool access easy.

- Meets or exceeds Americans with Disabilities Act Access Guidelines (ADA/ABA 2004)
- 300 pound lifting capacity
- Semi-permanent installation – easily removable from anchor
- Powered by a 24-volt rechargeable battery
- User operable via waterproof hand control
- 360° rotation ensures a safe transfer area
- Made from powder coated stainless steel and aluminum

Model No.	Description	Shipping Wt.	L	W	H
AXS1000	aXs	225 lbs.	40"	48"	48"

Includes battery, charger, battery console cover, waterproof control, stainless steel anchor socket with cover and spanner key, footrest and seat belt assembly

## Accessories:

- AX2320 Arm Rest Assembly
- AX9000 Total Cover
- 500-1000 Spineboard Attachment
- 900-2000 Stability Vest
- 100-1400 Wireless Controls

