



Arc Lite Rollator

Item #	Color	UOM
102EXL-BL	Blue	1/ea
102EXL-RD	Red	1/ea

- Innovative “Roman arch” frame design is lightweight while providing excellent support, stability and style.
- Center-folding cushioned seat is durable and waterproof.
- “X” style folding mechanism for compact storage & transport.
- Collapsible basket allows walker to be folded with the basket in place.
- Comfortable, form-fitting rubber handgrips.
- Durable, solid grip 8" caster wheels, ideal for both indoor and outdoor use.
- Low maintenance brake cables.
- Featherweight tires and wheels.

specifications

Handle Height	30.5" - 35.5"
Unit Dimensions	22.5" (W) x 26" (L)
Seat Height	21"
Seat Width	18"
Weight	15 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime