

BCW ADVANTAGE (850 LB. CAP.)



made in the USA

BCW ADVANTAGE

The **BCW Advantage** manual wheelchair is designed to meet the needs of the bariatric population. Constructed with a durable heavy-duty steel frame, four point telescoping seat and signature double x-brace support. The Advantage will fold for compact storage and has solid rubber tires and casters designed and tested for heavy-duty use. The upholstery is made with breathable fabric for added patient comfort. This wheelchair is available in a variety of sizes and will ship quickly to accommodate your needs.

FEATURES

- Adjustable seat-to-floor height
- Height adjustable and detachable desk-length armrest
- Quick ship
- 850 lb. weight capacity



Five year warranty on frame and x-brace.
Six months all other components
(excluding swingaway footrests).

AVAILABLE OPTIONS

Anti-tips • Seat belt • Telescoping IV pole • Oxygen tank holder
Elevating legrests • XL calf pads • LG footplates • **More options available**

OVERALL (DEPTH WITHOUT FRONT RIGGING)

SPECS

MODEL NO.	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	BACK HEIGHT	WIDTH	DEPTH	HEIGHT
AD24	24"	20"	19" or 20.5"	16"	31"	32"	40.5"
AD26	26"	20"	19" or 20.5"	16"	33"	32"	40.5"
AD28	28"	20"	19" or 20.5"	16"	35"	32"	40.5"
AD30	30"	20"	20.5" or 22"	16"	37"	32"	42.5"
AD32	32"	20"	20.5" or 22"	16"	39"	32"	42.5"
AD34	34"	20"	20.5" or 22"	16"	43"	33"	41"-42.5"

1-800-537-6454 www.wheelchairsofkansas.com

(For frame and upholstery color options see page 31)