

Accessories

MADE FOR COMFORT

Combo Cushion

Medline's combination seat and back cushion features high resiliency convoluted molded foam that incorporates Progressive Support Technology (PST) to ensure the proper degree of cushion firmness for the highest possible pressure relief. Hook and loop strap to secure cushion to wheelchair. Cover is removable, washable and flame retardant.

SIDEMINDER	SEAT W X D INCHES (CM)
MSCCOMB1616	16 x 16 (41 x 41)
MSCCOMB1816	18 x 16 (46 x 41)
MSCCOMB2018	20 x 18 (51 x 46)
MSCCOMB2218	22 x 18 (56 x 46)
MSCCOMB2418	24 x 18 (61 x 46)



Comfort Lap Tray

The Comfort Half Lap Tray is available as either an interlock or slide-on mount. The flip-up half lap pivots upward to allow users to easily enter and exit the wheelchair. With the slide-on mount option, the lap tray can be easily installed or removed. Also available in a full lap tray



TRAY	DESCRIPTION
MSCMLAPTRYL	Left Half
MSCMLAPTRYR	Right Half
MSCMEDLAPTRAY	Full Lap

SideMinder

As a fully adjustable lateral support, the SideMinder is recommended for users with moderate to severe impairment in their ability to maintain proper posture. This support is available with either standard or contoured lateral pads that adjust horizontally or vertically. For even the most difficult to correct postural asymmetries, the SideMinder is fully functional support.

SIDEMINDER	DESCRIPTION INCHES (CM)
MDTSDMDR	Fits 16-24 (41-61) wheelchairs



The Molded Foot

The Molded Foot offers a comfortable foot support designed with high quality polyurethane foam. It includes a rigid insert molded into the base for additional support. The molded foot can be attached to the wheelchair with straps or by bolting it directly to the foot plate.

MOLDED FOOT	DESCRIPTION
MSCMLDFT	Molded Foot



Medline Back Cushion

The Medline Back Cushion is a comfortable support cushion that fits in front of a standard wheelchair back, eliminating the hammock effect of a sling back. Two layers of high density foam provide support for the entire back. The ShearX all-way stretch cover reduces skin shear and is fluid-proof. Straps onto the wheelchair and is easy to remove.

HCPCS Code: E2611

BACK CUSHION	FITS WHEELCHAIR SIZE - INCHES (CM)
With Lumbar Support	
MSCBL16	16 (41)
MSCBL18	18 (46)
MSCBL20	20 (51)
Without Lumbar Support	
MSCB16	16 (41)
MSCB18	18 (46)
MSCB20	20 (51)



The Molded Arm

An ergonomically designed, cost-effective solution for users with arm issues. In addition to stabilizing the elbow, it features a convex palm area for user comfort as well as a gentle rise in elevation from elbow to fingertips help to prevent edema. Lightweight interlock bracket system locks down to either 3/4" or 7/8" tubing using a 2-bolt connection.

MOLDED ARM	DESCRIPTION
MSCMLDARML	Left Side
MSCMLDARMR	Right Side

