

Convoid Carrot 3 Child Restraint System

Canadian Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Carrot 3 Standard Features:

- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- Lower Anchors Connector System
- 5-Point Positioning Harness
- Trunk Positioning Inserts (Body Balance Pads)
- Seat Positioning Inserts
- Hip Positioning Inserts
- Chest Support: 9" Width
- Tether Strap (Black)

Standard Models:

<input type="checkbox"/> Carrot 3 CT3 BLACK - CCT3000-100	Price
<input type="checkbox"/> Carrot 3 CT3 BLUE - CCT3000-200	\$3,368
<input type="checkbox"/> Carrot 3 CT3 RED - CCT3000-300	\$3,368

All prices in CAD

Carrot 3 Recommended Item

<input type="checkbox"/> Top Securement Strap*	Price
<input type="checkbox"/> Lower Securement Strap*	\$79

*For vehicles not equipped with anchor points.

Trunk Positioning

<input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8")*	Price
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8")*	\$103
<input type="checkbox"/> Shoulder Protector Wings (Extends back height 6")	\$137
<input type="checkbox"/> Firm Lateral Trunk Supports - Black***	\$795
	\$282

*The 2" & 4" Backrest Extensions cannot be used simultaneously.

^Width adjustment from 10" - 12", and 2" of height adjustment.

**Manufacturer installation required.

Pelvic Positioning

<input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")]	Price
<input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")]	\$313
<input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")]	\$255
<input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")]	\$141
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	\$171
<input type="checkbox"/> Pommel - Black* (Medial Thigh Separator)	\$171
<input type="checkbox"/> Anti-Adduction/Knee Separator Pad	\$179
<input type="checkbox"/> Pelvis Lock Harness (Small) (20.5" - 33.75")	\$227
<input type="checkbox"/> Pelvis Lock Harness (Medium) (23" - 39.5")	\$377
<input type="checkbox"/> Pelvis Lock Harness (Large) (25" - 46")	\$377

^Growth specifications with listed extensions:

Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8".

*Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.

v7" Rodded Seat Extension required.

Foot & Leg Positioning

<input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")]	Price
<input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")]	\$484
<input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")]	\$484

*Not available with Rodded Seat Extension.

Small Support Kit

<input type="checkbox"/> Small Support Kit	Price
• Small Head Support: 7.25" x 8.5" (H x L)	\$694
• Small Chest Support: Width 6.9"	N/C

Accessories

<input type="checkbox"/> Upper Extremity Support (Tray) (Crash Tested)	Price
<input type="checkbox"/> 3D Body Balance Pads L/R Set	\$261
<input type="checkbox"/> Air-Mesh Cushions/Pads Set	\$227
<input type="checkbox"/> Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75")	\$415
<input type="checkbox"/> Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12")	\$241
<input type="checkbox"/> Anti-Escape Vest (Large) (Circum: 37.5" Height: 13")	\$241
<input type="checkbox"/> Removeable Incontinence Cover	\$241
<input type="checkbox"/> Hinge Lock Belt	\$302
<input type="checkbox"/> FAA Fixing Strap to Aircraft Seat	\$64
<input type="checkbox"/> Booster Kit	\$302
	\$1,736

Standard Shipping Rate: \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/27/2025

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0101-ENG-REV.00

Creating Possibilities

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"


^aBase model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Note: The same type of 2" & 4" extensions cannot be used simultaneously.

Inches
7"
2"

4"
2"

2"
4"
6"



Seat Depth Adjustments

2" Rodded Seat Depth Extensions

2" Rodded + 2" Seat Depth Extensions

7" Rodded Seat Depth Extensions

7" Rodded + 2" Seat Depth Extensions

7" Rodded + 4" Seat Depth Extensions

Position 1 (inch)	9.8"	11.8"	13.8"	16.8"	18.8	20.8"
Position 2 (inch)	11.3"	13.3"	15.3"	18.3"	20.3	22.3"

NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166