

Convoid EZ Rider - HCPCS E1236

Retail Price List / Order Form

Account #: _____

Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____

Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

EZ Rider Standard Features:

- Silver Powder Coated Frame
- Removable, Height Adjustable Swing-Away Footplates
- Tool-less Adjustable Seat Depth Growth of 5"
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (EZ12/14/16)
- 8" Front, 11" Rear Solid Tire (EZ18)
- Quick Release Rolko Rear Wheels (EZ12/14/16)
- One Piece Height Adjustable Push Handle (EZ12/14/16)
- Foot Operated Wheel Lock (EZ12/14/16)
- Zippered Pockets in Seat & Back Upholstery (Standard Upholstery Only: EZ12/14/16)

EZ Rider Convertible Features:

- Toggle Wheel Locks (EZ18)
- Two Piece Push Handles (EZ18)
- 2-Point Positioning Belt
- 10° Fixed Tilt
- Self-Tensio®
- Tool-less Upholstery Removal
- Tension Adjustable Back
- Made in the USA

Standard Models:	Price	Convertible Self-Propel Model	Price
☐ EZ Rider EZ12 - C900860	\$1,941	☐ EZ Rider EZ12 - C904469	\$2,566
☐ EZ Rider EZ14 - C900301	\$1,941	☐ EZ Rider EZ14 - C904470	\$2,566
☐ EZ Rider EZ16 - C900996	\$1,941	☐ EZ Rider EZ16 - C904471	\$2,566
☐ EZ Rider EZ18 - C900351	\$1,941	☐ EZ Rider EZ18 - C904472	\$2,566

- Upholstery Colors**
- | | Price |
|--|--------------|
| ☐ Standard Upholstery | |
| ☐ Panther Black | \$0 |
| ☐ Royal Blue | \$0 |
| ☐ Sassy Purple | \$0 |
| ☐ Princess Pink | \$0 |
| ☐ Forest Green | \$0 |
| ☐ Cherry Red | \$0 |
| ☐ Textilene Upholstery | |
| ☐ Navy Blue | \$0 |
| ☐ Violet Purple | \$0 |
| ☐ Wild Raspberry | \$0 |
| ☐ Teal Green | \$0 |
| ☐ Heavy Duty Reinforced Upholstery - Standard Upholstery [^] | \$280 |
| ☐ Heavy Duty Reinforced Upholstery - Textilene [^] | \$280 |
- [^]**Not available with Planar Option.**
***Heavy Duty Standard Upholstery is only available in Royal Blue, Sassy Purple, and Panther Black**

EZ12	EZ14	EZ16	EZ18
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855
C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465
C903467	C903467	C903467	C903467
C903466	C903466	C903466	C903466
C900835	C900321	C901031	C900392
C904494	C904495	C904496	C904497

- Frame Option**
- | | Price | HCPCS |
|--|--------------|--------------|
| ☐ Height Adjustable Flip-Up Armrests (pair)* | \$212 | E0973 |
| ☐ Two Piece Push Handle ¹ | \$0 | |
| ☐ One Piece Push Handle on Size 18 (Default is 2 piece) | \$900 | |
| ☐ Caster Lock For 2" Wide Front Wheels ¹ | \$98 | K0073 |
| ☐ Rear Anti-Tip Tubes (pair) | \$109 | E0971 |
- ¹**Manufacturer Installation Required.**
***See EZ Rider Convertible Options for item compatible with the Convertible Self-Propel Package.**

EZ12	EZ14	EZ16	EZ18
C904523	C904523	C904524	C904524
C904510	C904511	C904512	C904513
C902920	C902920	C902920	C902920
C900377	C900377	C900378	C900378

- Transit Required Items**
- | | Price | HCPCS |
|--|--------------|--------------|
| ☐ Transit Option ¹ | \$209 | |
| Head Support - Select Style Below: | | |
| ☐ Headrest Extension - Rectangular - Textilene** | \$135 | E0966 |
| ☐ Headrest Extension - Curved - Standard Upholstery (Not available in Planar) | \$135 | E0966 |
| ☐ Contoured Firm Headrest (Planar Option Only/Hrdw Req) | \$254 | E0955 |
| ☐ Adj. Removable Hardware (Req for Contoured Headrest) | \$257 | E1028 |
| Trunk Positioning - Select Style Below: | | |
| ☐ 5-Point Harness with Padded Covers (EZ12 only)* ¹ | \$150 | E0960 |
| ☐ H-Harness with Padded Covers (3-Pt Positioning Belt Required) | \$102 | E0960 |
| ☐ 3-Pt Positioning Belt | \$51 | E0978 |
| ☐ Full Torso Support Vest - Standard Upholstery (Not available in Planar) | \$275 | K0108 |
| ☐ Full Torso Support Vest - Textilene | \$275 | K0108 |
| Foot Positioners - Select Size Below: | | |
| ☐ Foot Positioners - Small 7" - 10" (pair) | \$163 | |
| ☐ Foot Positioners - Medium 8" - 11" (pair) | \$163 | |
| ☐ Foot Positioners - Large 9" - 12" (pair) | \$163 | |
- ¹**Manufacturer Installation Required.**
***Crash Tested - WC19 Compliant. Not available with Reinforced Upholstery, Cushion, Planar Option or H-Harness and 3 Pt. Belt.**
****Required when ordering the Contoured firm headrest and flat firm backrest extension.**

EZ12	EZ14	EZ16	EZ18
C900861	C900535	C900997	C900536
C900850	C900685	C901029	C900390
C902863	C900484	C900138	C902619
C904032	C904033	C904034	C904034
C904035	C904035	C904035	C904035
C903702			
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C900848	C900338	C901020	C905501
C905437	C905438	C905439	C900387
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

- Transit Options**
- | | Price | HCPCS |
|---|--------------|--------------|
| ☐ Q'Straint Transit Lap Belt | \$94 | E0978 |
| ☐ Sure-Lok Transit Lap Belt | \$94 | E0978 |

EZ12	EZ14	EZ16	EZ18
C903094	C903094	C903094	C903094
C903232	C903232	C903232	C903232

<u>Planar Option*</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ Planar Seat Modification ¹		\$359		C904049	C904050	C904051	C904052
☐ Planar Back with Attachment (Incl. 1 Backrest Ext. w/ Hrdw.) ¹		\$333	E2291	C904053	C904054	C904055	C904056
☐ Standard Upholstery with Planar Cushions (Default is Textilene)		\$440		C904900	C904900	C904900	C904900
Cushion - select one cushion below:							
☐ Support - General Use Cushion (Solid Seat Insert Included) ¹		\$251	E2601	C903684	C903687	C903690	C903693
☐ Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) ¹		\$486	E2607	C903686	C903688	C903692	C903695
☐ Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) ¹		\$486	E2607	C903685	C903689	C903691	C903694
Planar Head Support Options - Must Select One							
☐ Contoured Firm Headrest (Hardware Required) ¹		\$254	E0955	C904032	C904033	C904034	C904034
☐ Adj. Removable Hardware (Req for Contoured Headrest) ¹		\$257	E1028	C904035	C904035	C904035	C904035
☐ Flat Firm Backrest Extension (Hardware Required) ¹		\$238	E0955	C903928	C903582	C902151	C903674
☐ Adj. Removable Hardware (Req for Flat Firm Backrest Ext.) ¹		\$233	E1028	C904057	C904057	C904057	C904057
Planar Hard Lateral Supports							
Planar Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) ¹							
☐ Small = 3"x 5"		\$142	E0956	C902965	C902965	C902965	C902965
☐ Medium = 3.75"x 5.5"		\$142	E0956	C902969	C902969	C902969	C902969
☐ Large = 4.75"x 6"		\$142	E0956	C902973	C902973	C902973	C902973
Planar Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) ¹							
☐ Small = 3"x 5"		\$142	E0956	C904462	C904462	C904462	C904462
☐ Medium = 3.75"x 5.5"		\$142	E0956	C904463	C904463	C904463	C904463
☐ Large = 4.75"x 6"		\$142	E0956	C904464	C904464	C904464	C904464
☐ Planar Right Adj. Removable Hardware (Req for Lateral Support) ¹		\$185	E1028	C904030	C904030	C904031	C904031
☐ Planar Left Adj. Removable Hardware (Req for Lateral Support) ¹		\$185	E1028	C904030	C904030	C904031	C904031

¹Manufacturer Installation Required.^{*}Planar Option must be ordered with Planar Seat, Planar Back, Cushion, Planar Head Support, and Textilene Upholstery.

<u>EZ Rider Convertible Options</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ 8" x 2" Solid Front Tires		\$162				C903080	C903080
☐ Height Adjustable Flip-Up Armrests (pair)		\$212	E0973	C904044	C904044	C904045	C904045
☐ Upper Extremity Support Surface (Clear Tray) [^]		\$116	E0950	C903393	C903394	C903395	C903396
☐ Adj. Removable Hardware (Required for Clear Tray)		\$237	E1028	C903361	C903361	C903361	C903361

^{*}Requires Convertible Self-Propel Package.[^]Flip-Up Armrests and Hardware required.

<u>Head Supports</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ Padded Headwings [^]		\$254	E0955	C900849	C900682	C900683	C900684
☐ Occi Headwings [^]		\$238	E0955	C905271	C905271	C905271	C905272
☐ Headrest Extension - Rectangular - Textilene [^]		\$135	E0966	C900850	C900685	C901029	C900390
☐ Headrest Extension - Curved - Standard Upholstery [^]		\$135	E0966	C902863	C900484	C900138	C902619

[^]Not available with Planar Option.

<u>Trunk Positioning</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ H-Harness with Padded Covers (3-Pt Positioning Belt Req)		\$102	E0960	C903988	C903988	C903988	C903988
☐ Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Standard Upholstery [^]		\$142	E0956	C900846	C900327	C901018	C900386
☐ Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Standard Upholstery [^]		\$142	E0956	C900847	C900650	C900651	C900652
☐ Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Textilene [^]		\$142	E0956	C902859	C900678	C900140	C902614
☐ Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Textilene [^]		\$142	E0956	C902860	C900679	C900204	C902615
☐ Full Torso Support Vest - Standard Upholstery [*]		\$275	K0108	C900848	C900338	C901020	C905501
☐ Full Torso Support Vest - Textilene [*]		\$275	K0108	C905437	C905438	C905439	C900387
☐ Reducer Seat Insert (Reduces Seat Width 2" & Depth 1") [^]		\$368		C904211	C904110	C904110	C904110

¹Manufacturer Installation Required.[^]Not available with Planar Option.^{*}Not available with Single or Double Flap Lateral Support.

<u>Pelvic Positioning</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ 3-Pt Positioning Belt		\$51	E0978	C904017	C904017	C904017	C904017
☐ Support - General Use Cushion (Solid Seat Insert Included) ^{1*}		\$251	E2601	C903684	C903687	C903690	C903693
☐ Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) ^{1*}		\$486	E2607	C903686	C903688	C903692	C903695
☐ Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) ^{1*}		\$486	E2607	C903685	C903689	C903691	C903694
☐ Incontinence Cushion Liner (Cushion req) ¹		\$75	E2619	C902324	C902324	C902324	C902324
☐ Custom Seat Depth Extension 3" ^{^1}		\$98	E1297	C904003	C904005	C904007	C904009

¹Manufacturer Installation Required.[^]Not available with Planar Option.^{*}Not available with Seat Depth Extension or Standard Upholstery.

<u>Foot & Leg Positioning</u>		Price	HCPCS	EZ12	EZ14	EZ16	EZ18
☐ Medial Thigh Support-Abductor		\$169	E0957	C900843	C900329	C901024	C900388
☐ Lateral Thigh Support-Adductor		\$169	E0956	C900844	C900330	C901026	C900389
☐ Calf Panel [*]		\$64	E0995	C904580	C904580	C904580	C904580
☐ Angle Adjustable Footplates (pair)		\$166	K0040	C903268	C903268	C903269	C903269
☐ Footplate Securement Strap		\$32		C903584	C903584	C903584	C903584
☐ Foot Positioners - Small 7" - 10" (pair)		\$163		C903415	C903415	C903415	C903415
☐ Foot Positioners - Medium 8" - 11" (pair)		\$163		C903416	C903416	C903416	C903416
☐ Foot Positioners - Large 9" - 12" (pair)		\$163		C903417	C903417	C903417	C903417

^{*}Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

Tire Options

- ☐ 12.5" Solid Knobby Rear Tires*
☐ 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tires*
☐ 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tires*
☐ 8" x 2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tires*
☐ 8" x 2" Solid Front Tires
 Convertible Self-Propel Package
☐ 20" Mag Solid Tires
☐ 22" Mag Solid Tires
☐ 24" Mag Solid Tires
☐ Quick Release Axles
☐ Axle Attachment Brackets

Price

\$208
\$311
\$311
\$311
\$162
\$0
\$0
\$0
\$0
\$0

HCPCS

EZ12	EZ14	EZ16	EZ18
C904117	C904117	C904117	
C904118	C904118	C904118	
			C903084
			C903083
			C903080
C902940B			
	C902946B		
		C902946B	C902946B
C902942B	C902942B	C902942B	C902942B
C900377B	C900377B	C900378B	C900378B

***Not available with the Convertible Self-Propel Package.**

***With this option, standard tires are 6" Front, 8" Rear Solid Tires (EZ12/14) or 8" Front, 8" Rear Solid Tires (EZ16/18).**

Accessories

- ☐ Basic Headrest Cover (Canopy)
☐ Headrest Cover (Canopy)
☐ Extended Headrest Cover - No Windows (Canopy)
☐ Extended Headrest Cover - with Windows (Canopy)
☐ Upper Extremity Support Surface (Clear Tray)*^
☐ Adj. Removable Hardware (Required for Clear Tray)*
☐ Upper Extremity Support Surface (Foam Tray) Small^
☐ Upper Extremity Support Surface (Foam Tray) Medium^
☐ Upper Extremity Support Surface (Foam Tray) Large^
☐ Wheelchair Travel Bag
☐ Travel Bag for Large Wheels*
☐ Utility Bag
☐ Rain Cover Size 1 (Canopy Req)
☐ Rain Cover Size 2 (Canopy Req) (EZ14 Extended Canopy Req)
☐ Rain Cover Size 3 (Extended Canopy Req)
☐ Mosquito Net Size 1 (Canopy req)
☐ Mosquito Net Size 2 (Canopy req) (EZ14 Extended Canopy Req)
☐ Mosquito Net Size 3 (Extended Canopy req)

Price

\$180
\$238
\$260
\$307
\$116
\$237
\$158
\$158
\$158
\$215
\$135
\$40
\$157
\$157
\$157
\$157
\$157
\$157
\$157

HCPCS

EZ12	EZ14	EZ16	EZ18
C905021	C905022	C905023	C905024
C900853	C901210	C901205	C901211
C904721	C904725	C904726	C904727
C904737	C904741	C904742	C904743
C904525	C904526	C904527	C904528
C903369	C903369	C903369	C903369
C903669			
	C903670	C903670	
			C903671
C902022	C902022	C902022	C902022
C902484	C902485	C902486	C902487
C904531	C904531	C904531	C904531
C904939	C904939		
	C904940	C904940	C904940
			C905283
C904655	C904655	C904655	
	C904656	C904656	C904656
			C905319

***Manufacturer Installation Required.**

***See EZ Rider Convertible Options for item compatible with the Convertible Self-Propel Package.**

***Flip-Up Armrests and Hardware required.**

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The EZ Rider transit option includes four manufacturer installed wheelchair transportation anchors. The EZ Rider transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are EZ12: 66 lbs./30kg.; EZ14: 100 lbs./45.5kg; EZ16 and EZ18: 170 lbs./77kg.

The Convaid EZ Rider was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is being, or can be, made.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 4/28/2023

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF1200-ENG-REV.07

 Creating Possibilities