

## PreserveTech 360° Swivel Bath Chair

HCPCS: E0240,E0240

**PreserveTech™**  
Stays Cleaner • Lasts Longer



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

Compact and customizable, the PreserveTech Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility or other disabilities. Featuring a 360° swivel and locking at every 90°, the chair pivots for increased comfort and reduces friction wounds. The seat can be adjusted from 15 to 22, and the padded armrests can be raised out of the way or lowered and locked into place.

- Compact bath chair designed to aid in preventing friction and shearing injuries featuring 360 degree swivel action and locking at every 90 degrees . Eliminate the need for dangerous tub transfers caused by stepping over and in or out of the tub.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 350lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Padded armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Pivoting armrests can be raised out of the way.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15" to 22".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Push button, snap together design requires no tools and assembles in just minutes.
- Antimicrobial\* treated seat, backrest and non-slip tips, helps to minimize the growth of bacterial odors, mold, and fungi on treated surfaces.
- \*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

PRODUCTS	Description
RTL12A001-GR	PreserveTech Swivel Bath Chair, 1EA,1/CV
SPECIFICATIONS	Description
Material	Frame - Aluminum Seat Board - Polypropylene Backrest - HDPE Cushioning - EVA
Weight Capacity	8.80 lb
Seat To Floor Height	16"(W) x 18"(L)
Color	Grey
Dimensions	18"(W) x 16"(L) x 15-21"(H)
Retail Box Dimensions	7.09 "(W) x 20.08" (L) x 16.14" (H)
Retail Box Weight	10.30lb
Warranty Limited	Limited Lifetime
Weight without battery	350 lb
Width inside back legs	18"