

Drive DualEase Gel Cushion

HCPCS: E2622 (Pending),E2622 (Pending)




FIGURE A




FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Provide an elevated seating experience with the Drive DualEase Gel Cushion expertly engineered for pressure relief, skin protection, and personalized comfort. Perfect for individuals with limited mobility or anyone who values superior seated comfort, the Drive DualEase Gel Cushion adapts to your needs with a breathable gel design and effortless portability.
- Dual-layer unique polymer gel structure features a different firmness level on each side - simply flip the cushion to customize your comfort and relieve pressure.
- Flexible gel structure reduces shear and friction, supporting those with limited seated tolerance or at risk for pressure ulcers.
- Open-cell design and cooling hyper-elastic gel promote airflow and prevent moisture, consistently providing a cool seating experience.
- PreserveTech antimicrobial* stretch cover is low-shear, vapor permeable, and machine washable for easy care.
- Drive PreserveTech products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces. *Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.
- Built-in carry handle makes transport easy, while the non-skid bottom ensures the cushion stays put during use.
- 300lb. Weight Capacity
- 2-Year Warranty

PRODUCTS	Description
14992	Polymer Gel Cushion, 18x16x2, 1/ea, 4/cv
SPECIFICATIONS	Description
Warranty	Two Years
Material	Polymer Gel
Weight Capacity	300 lb.
Cover Material	Two-Way Stretch Nylon
Dimensions	14991: 16" x 16" 14992: 18" x 16" 14993: 20" x 16"
HCPCS	E2622 (Pending)
Weight	14991: 5 lb. 14992: 5.7 lb. 14993: 6.4 lb.