

AeroFly Ultra-Lite Wheelchair

HCPCS:



- Ultralight magnesium alloy frame weighs just 13.3 lbs. (excluding rear wheels, 19 lbs. with rear wheels)
- Includes both attendant brakes and wheel locks for assisted or independent use
- Fixed, desk length arms with molded armrest pads for comfort and durability
- Lightweight spoke wheels with treaded flat-free tires
- Quick release rear axles
- Fold down backrest
- Seat integrated carry handles
- Removable seat and back padded overlay cushions
- European style swing away footrests with ergonomic release levers
- Ships standard with rear anti-tippers for added safety

PRODUCTS	Description
AF18FDABK-SF	AeroFly,18",FixedDeskArm,SF,Black/Red
AF18FDAGY-SF	AeroFly,18",FixedDeskArm,SF,Grey/Blue
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	13.3 lbs. (Frame only weight, excludes rear wheels) 19 lbs. (Total weight inc. rear wheels)
Product Dimensions	30.3 x 26.8 x 36.6
Material	Magnesium Alloy
Seat Dimensions	18" x 16"
Weight Capacity	250 lbs.
Seat To Floor Height	19"
Wheels	22"
Carton Shipping Dimensions	31.1 x 13.4 x 30.7
Casters	6"
Folded Dimensions	30.3" x 13" x 29.7"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
AF18FDAGY-SF	AeroFly,18",FixedDeskArm,SF,Grey/Blue, EA