

## Universal Aluminum Crutch with Accessories

HCPCS: E0114



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Universal adjustment for use by a child, adult or tall adult (Figure A)
- Euro-Style clips and pushpin adjustment make it easy-to-adjust height of hand grip (Figure B)
- Lighter than wood, the aluminum crutch provides stability and durability
- Crutch disassembles for convenient storage (Figure C)
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Comfortable, durable underarm pad and hand grip

### PRODUCTS

### Description

|          |  |
|----------|--|
| RTL10433 | Universal Quick Adjust Crutch,1PR,8/CS |
|----------|--|

### SPECIFICATIONS

### Description

|                 |                  |
|-----------------|------------------|
| Warranty        | Limited Lifetime |
| Weight Capacity | 300 lb           |
| Height          | 38.5" - 62.5"    |