

Tilting Aluminum Rehab Shower Chair Commode

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- 30 Degrees of Tilting for Enhanced Comfort and Positioning in the Shower
- Versatile 3-in-1 design enables this model to easily be positioned over a standard toilet, used as a portable self-contained commode or shower chair.
- Simply and quickly adjust the height to achieve the ideal seated position.
- Padded flip-back armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Easy to empty commode bucket to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- Durable and easy-to-clean padded seat ,
- Adjustable Back Upholstery for Added Comfort
- Seat opening allows for unobstructed access to the perineal area for easy and effective cleansing.
- Swing-away detachable footrests with height and angle adjustments.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 300lbs
- Chest and Waist Strap Provide Safety and Comfort.

PRODUCTS

Description

NRS18010	Tilting Shower Chair Commode
----------	------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	17"(D) x 20"(W)
Weight Capacity	300 lb
Carton Shipping Weight	46lb
Weight	45 lb
Back Height	19"