

Quick 'N Easy Comfort Mattress

HCPCS: E0272-1,E0272-1



FIGURE A




FIGURE B



FIGURE C

- Durable high quality latex free polyurethane foam mattress for utmost comfort and support
- Non-slip bottom and 4 corner straps keep the mattress securely in place
- Two side carry handles make lifting and placement easier
- Compression technology creates a 14" x 14" x 37.5" box footprint creating convenience for transport and delivery
- PVC Polyurethane coated cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15076	Quick N Easy Comfort Mattress,1EA
15076C	Cover,15076,1ea
SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	275 lbs
FireRating	CFR 16 Part 1633 Federal Fire Code
Cover Material	PVC Polyurethane Coated
Dimensions	36" (W) x 80" (L) x 6" (H)
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]