

Compact Folding Rollator

HCPCS: E0143 + E0156,E0143 + E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

Constructed of durable aluminum, the Compact Folding Rollator is jam-packed with features to keep it a step ahead of the competition, including a 300 lbs weight capacity, small footprint for optimal shelf and storage space, and retail packaging to help your customers self-educate on product features and benefits. Introduce your customers to the perfect solution for maintaining their independence and mobility.

- Easy side-to-side fold for daily use and compact fold for storage and travel
- This lightweight rollator only weighs 15.4 lbs for effortless travel and transport
- Strong, durable aluminum frame supports up to 300 lbs
- Adjustable handle height from 31"-36" accommodates most users for a customized fit
- Push down on the handbrake to engage parking brake feature
- Padded backrest provides additional support when in a seated position
- Seat allows you to take a break when needed and remain comfortable
- 8" wheels are ideal for maneuvering both indoors and outdoors
- Removeable storage bag for your personal items
- Cane holder included for added convenience
- Tool-free assembly

PRODUCTS	Description
RTL10288BK	Compact Folding Rollator, Black, 1/EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	15.4 lb
Material	Aluminum
Seat Dimensions	17.95" x 8.66"
Handle Height	31" - 36"
Weight Capacity	300 lb
Seat To Floor Height	21.06"
Casters	8"
Color	Metallic Gray/Black
Dimensions	26" x 26.8"
Retail Box Dimensions	9.65" x 15.4" x25"