

Universal Seat Height, Aluminum Rollator with Fold Up and Removable Back Support, Padded Seat, 6" Casters with Loop Locks

Available
1 Pack


RTL10261BL

Blue, (Has replaced 301PSBN), 1/cs

RTL10261RD

Red, (Has replaced 301PSRN), 1/cs

NEW!



Blue Red

EASY HEIGHT ADJUSTMENT WITH SELF-THREADING KNOB!



10261RD-1

CASTER JOURNALS ARE OFFSET TO PROVIDE STRENGTH & DURABILITY!



2 ROLLATORS IN 1!
Seat Height Ranges from 18" to 22".



SEAT HEIGHT CAN BE ADJUSTED FROM 18" TO 22"!



FEATURES:

- Universal height adjustment on frame allows seat height to be adjusted from 18" to 22" in 1" increments.
- Easy height adjustment with self-threading knob.
- Comes with new seamless padded seat with zippered pouch under seat for added privacy and security of personal items.
- Removable, hinged padded back rest can be folded up or down.
- Caster journals are offset to provide strength and durability.
- Easy to use deluxe loop locks.
- Serrated brakes.
- Ergonomic handles are easy to grip and are adjustable in height.
- Weight capacity 300 lbs.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
RTL10261 SERIES	38"	29.5"	21"	25"	12"	13"	18" - 22"	14 lbs.	18 lbs.	300 lbs

drive™ rollators

To order call toll free: 877.224.0946