





Duogel




HIGH RISK

Features & Benefits

- Unique combination of TruGel and FloGel work together effectively to provide excellent pressure relief and weight distribution
- Gel and Foam combination allows for greater circulation and weight distribution
- Provides for medial and lateral thigh positioning as well as support for the greater trochanter of the thigh bone
- The thigh troughs are lined with TruGel and the Ischial area is lined with FloGel which is set on a viscoelastic foam giving greater comfort and seating stability
- Clear protective polyurethane film coating is easy to wipe clean and is suitable for incontinent users
- Lightweight and easy to transport
- Highly durable
- Includes outer cover



Pressure diagram with RELAX Duogel cushion



Pressure diagram without wheelchair cushion

Product Code	Dimensions	Weight
RD003-1514	38cm x 36cm x 8cm/ 15" x 14" x 3"	1.3kg/2.9lb
RD003-1616	40cm x 40cm x 8cm/ 16" x 16" x 3"	1.8kg/4.0lb
RD003-1618	40cm x 45cm x 8cm/ 16" x 18" x 3"	1.9kg/4.2lb
RD003-1717	43cm x 43cm x 8cm/ 17" x 17" x 3"	1.9kg/4.2lb
RD003-1816	45cm x 40cm x 8cm/ 18" x 16" x 3"	1.9kg/4.2lb
RD003-1818	45cm x 45cm x 8cm/ 18" x 18" x 3"	2.0kg/4.4lb
RD003-1820	45cm x 50cm x 8cm/ 18" x 20" x 3"	2.1kg/4.6lb
RD003-2018	50cm x 45 cm x 8cm/ 20" x 18" x 3"	2.1kg/4.6lb
RD003-2020	50cm x 50cm x 8cm/ 20" x 20" x 3"	2.5kg/5.5lb
RD003-2218	55cm x 45cm x 8cm/ 22" x 18" x 3"	2.8kg/6.2lb
RD003-2418	60cm x 45cm x 8cm/ 24" x 18" x 3"	3kg/6.6lb