

# 2025 order form

Use dimension chart to select appropriate size.


Scan for sizing and pricing:  
<https://www.rifton.com/e-pacer>



			K660		
* Frame includes one battery, one battery charger and charger power cord. (International power cords available on request.)			K661	\$5,480	<input type="checkbox"/>
* Color			CHOOSE	Red	<input type="checkbox"/>
				Blue	<input type="checkbox"/>
				Gray	<input type="checkbox"/>
				Lime	<input type="checkbox"/>
				Pink	<input type="checkbox"/>
Scale (must be chosen at the time of the original purchase. It cannot be added later.)			K662	\$1,300	<input type="checkbox"/>
Odometer			K663	\$220	<input type="checkbox"/>
Front handle			K665	\$415	<input type="checkbox"/>
Switch pole			K666	\$275	<input type="checkbox"/>
Arm supports (pair)	CHOOSE	Small arm prompts	K305	\$580	<input type="checkbox"/>
		Large arm prompts	K315	\$660	<input type="checkbox"/>
		Arm platforms	K397	\$435	<input type="checkbox"/>
		Arm platforms with handgrips	K396	\$490	<input type="checkbox"/>
Hip positioner	CHOOSE	Small	K667	\$360	<input type="checkbox"/>
		Large	K668	\$400	<input type="checkbox"/>
Hip positioner pad	CHOOSE	Small	K541	\$120	<input type="checkbox"/>
		Large	K542	\$160	<input type="checkbox"/>
Pelvic support	CHOOSE	Small	K306	\$295	<input type="checkbox"/>
		Medium	K316	\$330	<input type="checkbox"/>
		Large	K326	\$350	<input type="checkbox"/>
Thigh straps	CHOOSE	Narrow	K308	\$190	<input type="checkbox"/>
		Wide	K338	\$190	<input type="checkbox"/>
Ankle prompts			K684	\$240	<input type="checkbox"/>
Additional battery			K307	\$175	<input type="checkbox"/>
Additional thigh straps	CHOOSE	Narrow single	K318	\$95	<input type="checkbox"/>
		Narrow pair	K308	\$190	<input type="checkbox"/>
		Wide single	K328	\$95	<input type="checkbox"/>
		Wide pair	K338	\$190	<input type="checkbox"/>
Rifton accessories tote			K522	\$160	<input type="checkbox"/>

\* Required

# E-Pacer dimensions

User dimensions (inches)		K660
Girth		22 min–55 max
 <b>Key user dimension: girth</b>  <b>Important:</b> User's weight or height must not exceed: 350 lb or 77 inches		
Item dimensions (inches)		K660 E-Pacer
Overall length		44¾
Overall width		31½–51¾
Overall height		32¾–59
E-Pacer frame weight (lb)		72
Turning diameter		50
Floor to top of base leg for furniture clearance		8¾
Min. user armpit height (when standing or ambulating)		31½
Max. user height (when standing or ambulating)		77
Max. working load (lb)		350

## Pelvic support

User dimensions (inches)

Indicated by the red arrow below

Key user dimension:

half of thigh circumference

measured at the groin.

K306 small


K316 med.

K326 large

6-10

8-14

10-18



Product dimensions (inches)

K306 small

K316 med.

K326 large

a

8½

12

13½

b

16

20

25


c

9

12

14

## Hip positioner

Product dimensions (inches)–without pad		
		
	K667 small	K668 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11

The E-Pacer comes in your choice of five colors.



Red



Blue



Gray



Lime



Pink