

ENIGMA SPIRIT WHEELCHAIR



*Enigma Wheelchair
in Red shown,
item # XSES18RD*



**THE ULTIMATE IN
STYLE AND DESIGN!**



Enigma Spirit Wheelchair Single Axle

XSES18BL
XSES18RD

Blue, 1/cs
Red, 1/cs

FEATURES:

- Back folds down for storage and transport.
- Padded seat cushion.
- Adjustable tension back and seat.
- Puncture-proof solid front and rear tires.
- One touch adjustable armrests for easy transfer and comfort.
- Height adjustable anti-tip wheels.
- Palm activated quick release 24" rear wheels with black solid tires.
- Comes with swing-away removable footrests.
- New style ergonomic handgrips.
- Long reach brakes.
- Available in blue or red.
- Weight Capacity 250 lbs.
- Limited Lifetime Warranty



Red



Blue

SPECIFICATIONS

	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (w/riggings)	Product Weight (wo/riggings)	Carton Shipping Weight	Weight Limit
18"	28"	13"	36"	20"	17"	18"	16"	38"	29 lbs.	34 lbs.	250 lbs.

enigma spirit wheelchair

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM