

	EW-19	EW-20
Warranty	3 Year Limited Warranty	3 Year Limited Warranty
Colors	Blue & Black, Red & Black	Blue & Black, Red & Black
Weight Capacity	300 lbs	300 lbs
Fully Assembled	Yes	Yes
3 or 4 Wheel	3 Wheels	3 Wheels
FWD or RWD	Front Wheel Drive	Front Wheel Drive
Max Passengers	1	1
Airline Approved	No	No
Seat Dimensions (LxWxH)	10" x10" x 8"	16" x 13.5" x 3.25"
Type of Seat	Executive comfort with backrest	Large swivel seat with back rest
Type of Battery	Deep Cycle Sealed Lead Acid	Deep Cycle Sealed Lead Acid
Overall length	55"	55"
Overall width	26"	26"
Overall height	38"	38"
Armrests	No	Yes
Scooter Weight	128 lbs	128 lbs
Shipping Method	LTL Freight	LTL Freight
Top speed (Up To)	15 mph	15 mph
Distance (Up To)	21 miles	21 miles
Light(s)	LED Headlight, Rear Lights, Turn Signals Brake Lights	LED Front, Rear Lights, turn signals and brake lights
Turning Radius	53"	53 "
Braking Type	Hand Brakes	Hand Brakes
Left Hand Control Option	Yes	Yes
Horn	Yes	Yes
Incline capability	Up to 10 Degrees	Up to 10 Degrees
Ground Clearance	7"	7"
Suspension	Front and Rear	Front and Rear
Type of Tires	Front tubeless, rear tubed	Front tubeless, rear tubed
Battery Charger	48 volt offboard	48 volt offboard
Speed Control	Low-Med-High	Low-Med-High
Basket	Yes back	Yes back
Forward/Reverse	Yes	Yes
Throttle Type	Twist	Twist
Alarm	No	No
Front Tire Size	11" x 3"	11" x 3"
Rear Tire Size	11.5" x 3"	11.5" x 3"
Cell Phone Holder	Yes	Yes
USB Port	Yes	Yes
Voltage	48 volts	48 volts
Motor Wattage	500 watt	500 watt
Motor Type	Brushless Front Hub Motor	Brushless Front Hub Motor
Storage	Rear Basket	Rear Basket
Rear View Mirror	No	No
Extras	Folding Front Tiller, Removable Seat	Folding Tiller