

# RUBBER THRESHOLD RAMP



WITH BEVELED SIDES



Simply place the EZ-ACCESS Rubber Threshold ramp in front of the doorsill to provide a smooth ground-to-sill transition. This slip resistant ramp is designed to be used on any surface and can easily be trimmed or notched for a custom fit.

Made of 100% recycled materials. Ideal for indoor and outdoor use and can easily be moved to another location if needed. Ramp does not interfere with door tracks, making it ideal for use with sliding glass doors.

Usable width is 36" (42" overall including bevels) with a maximum height of 1½".

SPECIFICATIONS	THR36 150-1
Made of recycled rubber	•
Non-skid surface	•
Usable size (H x L x W)	 1½" x 12" x 36"
Overall size (H x L x W)	 1½" x 12" x 42"
Weight capacity (lbs.)	850
For wheelchairs & scooters	•
Usage guidelines	1½" rise



Bariatric weight capacity = 800 lbs. or greater



Bariatric width = 36" or greater



 UPS Shippable

Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.