



LC-358M shown in  
Cloud 9 Stone fabric

MADE IN  
**USA**  
OF FOREIGN & DOMESTIC COMPONENTS

# HERITAGE COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS

Comfort and style complement one another perfectly in the **Heritage Collection** from Pride® Power Lift Recliners. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection delivers style and durability that will last for years to come. Pride's Heritage Collection offers numerous standard features, superb comfort and quiet and smooth performance, all at a great value.

(US) 800-800-8586 | (CAN) 888-570-1113  
[pridemobility.com](http://pridemobility.com) | [prideliftchairs.com](http://prideliftchairs.com)

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —



Available on the Heritage Collection Series (except 558, 358XL & XXL)



SPECIFICATIONS	Medium	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
	LC-558 <sup>8</sup>	LC-358PW <sup>2,5</sup>	LC-358S <sup>5</sup>	LC-358M <sup>2,5</sup>	LC-358L <sup>2,5</sup>	LC-358XL <sup>4,6</sup>	LC-358XXL <sup>4,6,7</sup>
	Infinite, Lay Flat Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'8" - 6'2"	3-Position Chaise Lounger 5'6" - 6'1"	2-Position Chaise Lounger 5'6" - 6'1"
Weight Capacity	325 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	19"	19"	19"	20"	21"	20"	20"
Seat Depth	20"	19"	19"	20"	22"	22"	22"
Seat Width	19.5"	24"	18"	20"	22"	26"	30"
Top of Back to Seat	26"	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage	NA	Optional	Optional	Optional	Optional	Optional	NA
Seating Type	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Sinuous Spring	Sinuous Spring
Distance From Wall	24"	18"	18"	17"	19"	19"	14"
Overall Chair Width	35"	40"	35"	35"	37"	39.5"	43"
Chair Weight	147 lbs.	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Warranty	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated	7-Year Prorated

## HERITAGE FABRIC SELECTION

### CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT

### LEXIS STA-KLEEN®

sta\*kleen.  
PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

### ULTRALEATHER™



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

### ARIA<sup>3</sup>



ESPRESSO



COOL GREY



RED



LAZULI



SAND

- Batteries not included.
- Standard footrest extension.
- Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
- Removable, fully padded chaise.

- Unique spring design and a one-piece soft foam seat for outstanding relaxation.

- Features Dual Motor Sync System.

- Only available in Stone, Walnut, Lexis Sta-Kleen, Chestnut and Mushroom.

- Standard with footrest extension and lumbar pillow.

- Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

- Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Look for the Quick Ship Indicator for the shortest lead times available.

NEW!

## FEATURES AND BENEFITS

- Standard footrest extension on 358PW, M and L for additional leg support while reclined
- Easily adjustable leg levelers provides additional support for uneven surfaces
- Easy-to-use hand control with USB charger port (not available when adding optional heat and massage)
- Standard head and arm covers
- Large dual pockets for extra storage
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-358XL and LC-358XXL)
- Quiet and smooth lift and scissor mechanism system
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Premium furniture grade laminate/hardwood frame for added strength and durability

## OPTIONS\*

- Heat and massage (except XXL and LC-558)
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Wireless remote (not available on LC558, XL or XXL)
- FirstUp™ technology (except LC-358XL, LC-358XXL, and LC-558)
- Non-skid leg levelers

Additional lead time may apply for some options. \*Additional charges apply for all options.

(US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —