

Hybrid 2 | SPECIALTY

- Combination wheelchair and transport chair
- Two chairs in one
- Rear wheels remove with just the push of a button
- Weighs only 33 lbs (15 kg), 24 lbs (11 kg) in transport
- Converts easily in seconds
- Comfortable nylon upholstery
- Includes anti-tippers

ACCESSORIES

MDS85181U

O₂ Holder

MDS85196

Anti-Fold Device

MDS85183

I.V. Pole

MDS85197A

Overhead Anti-Theft Device

WCA806991

Tinnerman Legrest Locks

MDS85190

O₂ / I.V. Combo



MDS806250H2

HYBRID 2 IN TRANSPORT MODE

CHAIR	ARM STYLE	LEG STYLE	SEAT W X D INCHES (CM)	HCPCS*	WT. CAP. LB (KG)
MDS806250NH2	Desk Length Removable	Swing-Away Detachable Footrest	16 x 16 (41 x 41)	Pending	300 (136)
MDS806150H2	Full Length Permanent	Swing-Away Detachable Footrest	18 x 16 (46 x 41)	Pending	300 (136)
MDS806200H2	Full Length Permanent	Swing-Away Detachable Elevating Legrest	18 x 16 (46 x 41)	Pending	300 (136)
MDS806250H2	Desk Length Removable	Swing-Away Detachable Footrest	18 x 16 (46 x 41)	Pending	300 (136)
MDS806300H2	Desk Length Removable	Swing-Away Detachable Elevating Legrest	18 x 16 (46 x 41)	Pending*	300 (136)

*Medline submitted application for K0004 code