

Part Number	Jazzy® Carbon HD
Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ³	1" at motor
Turning Radius ³	36.4"
Overall Length ³	36.81"
Overall Width ³	26.4"
Suspension	Full front and rear suspension
Tires	Front: 6.5" Rear: 10.6"
Range Per Charge ^{1,2}	Up to 12.4 miles per charge
Product Weight ^{3,5}	32.2 lbs. without batteries
Standard Seating ³	Dimensions: 20" x 20" seat Seat sling with removable cushion
Folded Dimensions	31.3" L x 26.4" W x 13" H
Motor Output	24V 250W Brushless
Brakes	Regenerative and electromechanical
Battery Requirements	Type: Lithium battery box Size: 20 Ah
Battery Weight ⁴	7.5 lbs
Battery Charger	2 Amp
Warranty	5-year on frame; 1-year limited on electronics 1-year limited on drivetrain 1-year limited on battery

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- Battery weight may vary based on battery manufacturer.
- Includes standard seat weight.

* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

Jazzy®

carbon HD

Weighs just
32lbs
400lbs
weight capacity

PRIDE MOBILITY®

800-800-8586 • pridemobility.com

FDA Class II
Medical Device

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.


081425



Folds and
unfolds for easy
transport



Quick disconnect
joystick with
LCD display



Full Front
and rear
suspension

Jazzy[®] carbon HD

Weighs just
32lbs
400lbs
weight capacity

Quick disconnect
LCD joystick with
digital display, USB-C
port & puddle light

Shown with optional
headrest and attendant
control kit

20" x 20" seat with
padded back and
storage pocket

Mesh cup holder
(Standard)

Adjustable lap belt
(Standard)

Compact base
folds & unfolds
for easy transport

Durable, lightweight
carbon fiber frame

Height adjustable
foot platform

Full front &
rear suspension

FDA Class II
Medical Device*



Folds & unfolds
for easy transport



Quick disconnect joystick
with LCD display



Full front &
rear suspension



Latch for locking
unit open or closed

HEAVY-DUTY. LIGHTWEIGHT.

The Jazzy[®] Carbon HD weighs 32.2 lbs. without the battery and boasts a 400-lb. weight capacity without compromising on the quality or performance. It quickly folds and is easy to pick up thanks to a durable, lightweight carbon fiber frame.

FEATURES

- Unit weighs 32.2 lbs. without the battery
- 400-lb. weight capacity
- Durable carbon fiber frame
- Front and rear suspension
- Top speed of up to 3.7 mph
- Folded dimensions: 31.3" L x 26.4" W x 13 H
- Easy fold design that locks into folded position
- Flip up lightweight carbon fiber arms
- Padded seat back with removable cushion
- Extra-wide adjustable foot platform 17.25" W x 5.9" D
- 20 Ah lithium battery** charges quickly (Standard)
- Storage pocket on back of seat
- Quick disconnect LCD joystick with digital display, USB-C port & puddle light
- Mesh cup holder, cushion and adjustable lap belt, (standard)

COLOR OPTIONS



OPTIONAL ACCESSORIES

Additional 12Ah
airline-compliant battery**

**Additional 20Ah
lithium battery (Coming Soon)**

Battery Travel Bag

Headrest with Attendant
Control Bracket

Jazzy[®] Travel Kit (Coming Soon)

Military Patches

Simplicity Cushion 20"x 20"

**Elevating Leg Rests
(Coming Soon)**



Optional Jazzy[®] Travel Kit
(Coming Soon)

**Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Consult your airline prior to flight to ensure your lithium battery is permitted.