



# Hannah 135 with Full electric Platform



Designed as a tall rail bed and built extremely stable, with a full electric platform, that offers electric height adjustment and electric head and knee elevation. Hannah has boundaries for the active patient and positioning when needed for care and medical issues.

*Available with options to meet every need.*

## Measurements:

- Safety rail height 53"
- Weight capacity 330 pounds
- 18.5" transfer height
- Extra long twin mattress
- Inside dimensions: 79" X 39.5"
- Outside dimensions: 85" X 44"
- Overall height: 72"

## Specifications:

- Solid natural beech wood
- Sealed for disinfecting
- Special door locking system
- Low transfer height
- Mattress sits above the bed frame
- 4 doors on one side
- Wood rods allow air circulation
- FDA Registered




REMOVABLE PADDING    NON-REMOVABLE PADDING



WOOD TOP

MESH FABRIC TOP



STANDARD MATTRESS

DURABLE COVER MATTRESS



PLEXIGLAS OVER RODS