



Hannah 135 with Manual Wood Platform



Designed as a tall rail bed and built extremely stable, with a wood platform that offers manual head and foot elevation. Hannah has boundaries for the active patient, who needs a safe environment to protect them from harm.

Available with options to meet every need.

Measurements:

- Safety rail height 53"
- Weight capacity 330 pounds
- 18.5" transfer height
- Extra long twin mattress
- Inside dimensions: 79" X 39.5"
- Outside dimensions: 85" X 44"
- Overall height: 72"

Specifications:

- Solid natural beech wood
- Sealed for disinfecting
- Special door locking system
- Low transfer height
- Mattress sits above the bed frame
- 4 doors on one side
- Wood rods allow air circulation
- FDA Registered




REMOVABLE PADDING NON-REMOVABLE PADDING



WOOD TOP

MESH FABRIC TOP



STANDARD MATTRESS

DURABLE COVER MATTRESS



PLEXIGLAS OVER RODS