



## Hannah 170 with Manual Wood Platform



Designed as a tall rail bed and built extremely stable, with a wood platform that offers manual head and foot elevation. Hannah has boundaries for the active patient, who needs a safe environment to protect them from harm.

*Available with options to meet every need.*

### Measurements:

- Safety rail height 67"
- Weight capacity 330 pounds
- 18.5" transfer height
- Extra long twin mattress
- Inside dimensions: 79" X 39.5"
- Outside dimensions: 85" X 44"
- Overall height: 86"

### Specifications:

- Solid natural beech wood
- Sealed for disinfecting
- Special door locking system
- Low transfer height
- Mattress sits above the bed frame
- 4 doors on one side
- Wood rods allow air circulation
- FDA Registered

