

AirLite Pro

Personal Mobility Vehicles from Küschall®



**Go ahead,
pinch yourself...**

Ultra lightweight design

Dialed in driver adjustments

High performance options

Titanium where titanium counts

Küschall
Champion



Innovative Design Features

- Center mass adjustment
- Fully adjustable seat module
- Folding back frame
- Titanium footrest, seat and back frame
- Carbon fiber camber tube
- 75° & 90° front frame angles
- Dual cylinder suspension system
- Vario-Ax adjustable camber tube
- High performance rear wheels
- STAREC aluminum front casters

Vehicle Specifications

Weight capacity: 267 lbs.

Vehicle weight: 19 lbs. as shown

Seat widths:

>13"/33cm to 18"/45cm

Seat depths:

14"/35cm to 20"/50cm

Toll-free: (888) 682-2571
www.kuschallna.com

küschall®
designed for life

Compact

Personal Mobility Vehicles from Küschall®



**It's more than
meets the eye...**

Amazingly ultra-light
Sleek innovative design
Fully adjustable
Folds compactly

Küschall
Ultra-Light
▶

Küschall
AirLite Pro
▶

Küschall
Champion
▶



Innovative Design Features

- Hemi height standard frame
- Hemi height abduction frame
- Fully adjustable axle bracket
- Posi-lock seat frame
- Adjustable seat upholstery
- Spring loaded wheel locks
- Anti-rotation armrest bracket
- Painted oval cross tubes
- Integrated caster flow system
- High performance rear wheels
- STAREC aluminum front casters
- 70° & 80° hangers swing in & out

Vehicle Specifications

Weight capacity: 289 lbs.
Vehicle weight: 24 lbs. as shown
Seat widths:
 <13"/32cm to 20"/50cm
Seat depths:
 <13"/32cm to 20"/50cm

Toll-free: (888) 682-2571
www.kuschallna.com

küschall®
designed for life

Ultra-Light

Personal Mobility Vehicles from Küschall®

**Think "outside"
the box...**

Unbelievably ultra-light

High performance feel

Smooth streamlined design

Adjusts for a personalized fit



Küschall
Compact
▲



Innovative Design Features

- Fully adjustable axle bracket
- Posi-lock seat frame
- Adjustable seat upholstery
- Spring loaded wheel locks
- Anti-rotation armrest bracket
- Painted oval cross tubes
- Integrated caster flow system
- High performance rear wheels
- STAREC aluminum front casters
- 75° or 90° fixed front frames
- Aero back upholstery

Vehicle Specifications

Weight capacity: 289 lbs.

Vehicle weight: 22 lbs. as shown

Seat widths:

<13"/32cm to 20"/50cm

Seat depths:

<13"/32cm to 20"/50cm

Toll-free: (888) 682-2571
www.kuschallna.com

küschall®
designed for life

Champion

Personal Mobility Vehicles from Küschall®

**Finally, the
decision is easy...**

**Rigid vehicle performance
Folding vehicle convenience**



Küschall
Compact ▲

Küschall
Ultra-Light ▲

Küschall
AirLife Pro ▲



Innovative Design Features

- One hand folding capability
- Horizontal scissor lock
- Auto-up footrest plate
- Fold down back rest
- Fully adjustable axle bracket
- Carbon fiber frame option
- 75° & 90° front frame angles
- High performance rear wheels
- STAREC aluminum front casters

Vehicle Specifications

Weight capacity: 267 lbs.
Vehicle weight: 24 lbs. as shown
Seat widths:
 >13.5"/34cm to >17.5"/44cm
Seat depths:
 >13.5"/34cm to <18.5"/46cm

© 2004, Küschall. 10.04.

küschall®
designed for life

Küschall
1831-A Lefthand Circle
Longmont, Colorado 80501

Toll-free: (888) 682-2571
Fax: (866) 651-6973
www.kuschallna.com