

INVACARE®

Leo™

Scooter Series



Yes, you can.®

INVACARE®

Leo™ Scooter Series



The Invacare® Leo™ 3-Wheel and 4-Wheel Midsize Scooters offer a smooth, safe drive that is easy to maneuver.

Features include a 350 lb. weight capacity, standard full lighting package, built-in splash guards to protect the electronics and transaxle, comfortable seating that swivels and slides and flat free tires. With a sophisticated design and impressive features, the Leo Scooter will **get you where you're going in style!**

Features



Quick Adjust Tiller

Infinitely adjustable tiller allows user to suit their needs. Charger port is conveniently located on the tiller.

Intuitive Control Panel

Offers large buttons that are easy to operate and bright lights that are easy to distinguish.

Full Lighting Package

Standard LED headlights and rear brake lights provide increased visibility.

Note: Do not operate on roads, streets or highways.



Accessories



Rear View Mirror - Model no. ACC250
Removable mirror can be mounted on the left or right.



Lockable Front Box - Model no. ACC260/ACC270
Secure box for valuables in blue or silver to match the scooter. Dimensions are 17 "L x 9 "W x 11.5 "H.



Rear Basket - Model no. ACC230
Large basket is ideal for storing extra items. Rear Basket dimensions are 18 "L x 13 "W x 8 "H.



Safety Flag - Model no. ACC140
Bright orange flag makes the scooter easy to see from a distance.



Crutch and Cane Holder - Model no. ACC200
Keeps a crutch or cane stable during transit and within reach.



Oxygen Holder - Model no. ACC220
Safe storage of oxygen during transit.



Walker Holder - Model no. ACC210
Keeps a walker secure during the drive.

Specifications

	Invacare® Leo™ 3-Wheel Scooter	Invacare® Leo™ 4-Wheel Scooter
Model No:	LEO-3B, LEO-3S	LEO-4B, LEO-4S
Overall dimensions		
Base length	47.2"	48"
Base width	23.2"	23.2"
Ground clearance	2.4"	2.5"
Seat dimensions		
Width	19"	19"
Depth	16"	16"
Back height	16.5" - 19.5"	16.5" - 19.5"
Arm dimensions		
Height	8"	8"
Speed	5 mph	5 mph
Weight limit	350 lb.	350 lb.
Maximum incline capability	6° grade	8° grade
Battery	U1	U1
Estimated travel range*	18.6 miles	15 miles
Tires	10"	10"
Turning radius	47.2"	52"
Weight	173.7 lb.	184 lb.
Accessory weight capacity		
Lockable front box	10 lb.	10 lb.
Rear basket	10 lb.	10 lb.
Limited warranty	1 year	1 year
Frame, electronics, transaxle		

*Travel range will vary depending on user weight, speed, battery charge and terrain.

Shroud Colors



Blue

Silver

Disassembles for transport.



Note: Do not operate on roads, streets or highways.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2011

Invacare Corporation.

All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.

Form No. 09-044
Rev. 07/11
110621

