



## Liko® M220 and M230 Mobile Lifts

**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™

# Designed for Residents, Caregivers, and Decision Makers

## Simply, easy to use

The easy-to-use Liko® M220 and M230 Mobile Lifts are designed to meet lifting needs in post-acute care environments with the durability and quality you've come to trust from Hill-Rom. Partnering with our safe lifting expertise can help your facility minimize financial risk from caregivers and resident injuries.

The Liko M220 and M230 Mobile Lifts feature intuitive design elements that make them easy to use while meeting the daily lifting needs in your facility. The lifts create a secure and comfortable lifting experience for both residents and caregivers.



The Liko M220 and M230 Mobile Lifts were designed to compliment the full portfolio of Hill-Rom® post-acute care products including bed frames, prevention and therapy surfaces, bariatric solutions and furniture.



### Residents

- Full electric lifting motion
  - High/low
- Gentle speed for patient safety and comfort
- Manages all daily lift situations
  - To/from chair
  - To/from bed
  - To/from toilet
  - To/from floor
  - To/from bath

### Caregivers

- Intuitive controls
  - Easy to use and understand
- Ideal for lifting needs in an extended care environment
- Electric lift mechanism allows caregiver to concentrate on patient (rather than pumping or cranking standard hydraulic/manual lifts)
- Long lift arm
  - Easy access to middle of bed
  - Allows comfortable rotation
  - Lifting height of 67.7 inches

### Decision-Makers

- Minimize financial risk
- Resident and caregiver safety
- Resident and caregiver satisfaction

# Liko® M220/M230 Mobile Lift Features

## Electric high/low lifting motion

- Provides safety and comfort to resident and caregiver

## Dual-mode emergency lowering

- Electrical emergency lowering uses available battery power to safely lower resident
- Manual emergency lowering features knob with arrows showing direction to turn

## Easy-to-use hand control

- Cable length allows caregiver to maintain contact with resident while using lift
- M230 four-button hand control for lifting and base width adjustment
- M220 two-button hand control for lifting

## Quick reference guide

- Reminds caregivers how to operate lift

## Long-lasting battery

- Charged while on the lift or easily removable for optional wall charging

## Emergency stop button

- Quickly and safely stops electrical functions

## Low battery indication

- Audible battery indication while pressing hand control button signals caregiver when battery is low

## Long Lift Arm

- Provides access to center of bed
- Facilitates daily resident transfers and lifting from floor

## Easy-to-use, standard two-point sling bar

- Lifts up to 400 lbs
- Compatible with Liko slings, scales and accessories

## Gentle actuator speed for patient comfort and safety

## Low-profile 4.1" leg clearance to fit easily under beds

- Front 3" twin wheels and rear 3" twin wheels with brakes are non-marring

## Liko M230 electric base

- Allows caregivers to access resident in a variety of lifting situations with the touch of a button

Liko M220 features a manual base with easy-to-use foot pedal for adjustments

## Technical Specifications

### Liko M230 (electric base)

Lifting capacity.....	400 lbs
Lifting interval.....	22.4" to 67.7"
Wheels.....	Front 3" twin wheels, Rear 3" twin wheels with brakes
Material.....	Steel
Emergency lowering.....	Manual and Electrical
Battery.....	In series: 2x12V 2.9 Ah valve-regulated lead-acid gel batteries
Battery charger.....	Built-in, 100-240V AC, 50-60 Hz, max 400 mA
Total unit weight.....	85 lbs
Heaviest removable part.....	44.75 lbs

Learn more at  
[www.hill-rom.com/m230](http://www.hill-rom.com/m230)

Includes additional product details  
 and safe lifting information.

# Liko offers a wide array of accessories for the M220/M230 Lift.



## Slings

For more than 30 years, we have been studying lifting situations in different settings. This knowledge has given us the most functional, safest and complete range of slings on the market.



## Liko® 400 Scale

A digital scale that can be used in combination with all of Liko's lifts. Connected above the slingbar, the Liko Scale is easy-to-use and read.



## Charger Station

With optional wall-mounted charger station and extra battery pack, reserve power is always available.

## We have the safe lifting solution for your long-term care needs.

Since 1929, Hill-Rom has committed itself to providing high-quality, innovative products and solutions, along with support from a dedicated team of clinical experts and service personnel.

Today, our Extended Care Team is proud to offer a full line of Liko® safe lifting solutions for caregivers and residents through our trained sales force and national network of Authorized Liko Distributors.

With more than 30 years of experience, Liko is a global leader in developing safe lifting solutions for many different care and work environments, including long-term care. Liko authorized distributors are trained in moving and handling patients and can partner with you to develop a safe lifting program for your system or facility.

You can count on the people, processes and technology of Hill-Rom—in hospitals, extended care facilities and in the home.



Mobile Lifts



Sit-to-Stand Lifts



Overhead Systems

Liko® M220/M230 has been tested by an accredited testing institute and complies with the requirements for Medical Device Directive Class I products (MDD 93/42/EEC). Liko® M220/M230 complies with the requirements in EN ISO 10535, IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, ANSI/AAMI ES60601-1, and CAN/CSA C22.2 no 60601-1. Design and Quality by Liko in Sweden. Liko® is quality certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485. Liko® is certified in accordance with environmental standard ISO 14001.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2011 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

**Hill-Rom**