## 2024 order form for medium

Use dimension chart to select appropriate size.



To order or get a quote: Go online: https://www.rif Email: sales@rifton.com	ton.com	/activity-sizing		L SALIT				
<b>Call:</b> 800.571.8198			Standard R840 medium		1	Hi/lo R850 medium		d
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	운	Adjustable	R841	\$670		R851	\$1390	
	OSE	Adjustable w/ back spring	R842	\$760		R852	\$1510	
* Base	웆	Without spring	R843	\$385		R853	\$2790	
	OOSE	With spring	R844	\$465		N/A	•	
* Arm supports (pair)	CHC	Armrests	R868	\$270		R868	\$270	
	3SO(	Forearm prompts	R869	\$675		R869	\$675	
* Adjustable legs		Short	R893	\$65		N/A		
	운	Long	R894	\$85		N/A	•	
	CHOOSE	Short with casters	R895	\$330		N/A	•	
		Long with casters	R896	\$360		N/A	•	
* Pads		Tan	R835	\$275		R835	\$275	
		Blue	R836	\$275		R836	\$275	
	용	Red	R837	\$275		R837	\$275	
	CHOOSE	Pink	R974	\$275		R974	\$275	
		Purple	R975	\$275		R975	\$275	
		Green	R976	\$275		R976	\$275	
Headrest	0	Contoured	R897	\$260		R897	\$260	
	CHOOSE	Flat	R898	\$260		R898	\$260	
	m	Adjustable winged	R881	\$290		R881	\$290	
Laterals (pair)		Small	R828	\$220		R828	\$220	
	CHOOSE	Small with chest strap	R829	\$275		R829	\$275	
	m	Large	R838	\$235		R838	\$235	
Push handles (recommended for all chairs with wheels)			R878	\$166		Included		
Tray (requires armrests)			R668	\$300		R668	\$300	
Handhold			R649	\$31		R649	\$31	
Chest strap (wide)			R863	\$100		R863	\$100	
Butterfly harness	운	Regular	R864	\$135		R864	\$135	
	OOSE	Slim cut	R882	\$135		R882	\$135	
Thigh belt			R865	\$105		R865	\$105	
Pelvic harness			R866	\$190		R866	\$190	
Hip guides (pair)			R867	\$230		R867	\$230	
Abductor	윤	Small	R884	\$140		R884	\$140	
	CHOOSE	Large	R891	\$160		R891	\$160	
Adductors (pair)			R892	\$220		R892	\$220	
Leg prompt			R858	\$245		R858	\$245	
Footboard (recommended for chairs with long legs)			R848	\$525		Included		
Ankle straps (pair)			R899	\$100		R899	\$100	
Sandals (pair)			R954	\$105		R954	\$105	
Wedges (pair, require sandals)			R698	\$80		R698	\$80	
Lumbar and seat support kit			R887	\$88		R887	\$88	
Backrest filler pad			R849	\$135		R849	\$135	
Whitmyer headrest adapter			R804	\$63		R804	\$63	



40-56

## **Room to grow**

## Medium chair at its smallest dimensions

- Backrest height 15½"
- 2 Distance between laterals 61/2"
- Seat depth 11"
- 4 Seat width 81/2"
- Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



# Medium chair at its largest dimensions

- Backrest height 19"
- Distance between laterals 12"
- Seat depth 16"
- 4 Seat width 14"
- ⑤ Armrest height above seat 10½"
- 6 Seat height above footboard 151/4"





### **Medium standard base**

R840 Rifton Activity Chair

#### **User dimensions (inches)**

Height Key user dimension: height



The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)	R840 medium Standard base
Frame width	short legs: 231/4
	long legs: 241/2
	short legs w/ casters: 241/2
	long legs w/ casters: 251/2
Seat height above floor	short legs: 12½-16½
	long legs: 19-23
	short legs w/ casters: 16-20
	long legs w/ casters: 19-23
Seat angle (tilt-in-space) – degrees	15° forward, 15° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	73/4-151/4
Seat width with hip guides without hip guides	8½-11½ 14
Seat depth	11-16
Armrest height above seat	7-101/2
Distance between laterals	6½-12
Backrest height	15½ – 19
Headrest height above seat	17½-24½
Max. working load (lb)	150



### **Medium hi/lo base** R850 Rifton Activity Chair

### User dimensions (inches)

Height



**Key user dimension: height**The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)	R850 medium Hi/lo base
Frame width	271/2
Seat height above floor	12-23
Seat angle (tilt-in-space) – degrees	15° forward, 25° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45°-110°
Seat height above footboard	73/4 – 151/4
Seat width with hip guides without hip guides	8½-11½ 14
Seat depth	11-16
Armrest height above seat	7-101/2
Distance between laterals	61/2 – 12
Backrest height	15½-19
Headrest height above seat	171/2-241/2
Max. working load (lb)	150

