

2026 Tricycle order form

Use dimension chart to select appropriate size.

		X320 small	X330 medium	X340 large
* Tricycle frame includes pedals, parking brake with variable resistance, and a rear basket.		X321 \$1,950 ☐	X331 \$2,230 ☐	X341 \$2,910 ☐
* Color (Note that the choices are different between sizes.)		CHOOSE Blue ☐ Red ☐ Teal ☐ Lime ☐ Pink ☐ Blue ☐ Red ☐ Teal ☐ Lime ☐ Pink ☐ Blue ☐ Red ☐ Teal ☐ Black ☐ White ☐		
* Handlebar	CHOOSE Loop	X322 \$235 ☐	X332 \$245 ☐	X342 \$245 ☐
	Conventional	X323 \$235 ☐	X333 \$245 ☐	X343 \$245 ☐
* Back support	CHOOSE Small backrest pad	X326 \$165 ☐	X326 \$165 ☐	N/A
	Large backrest pad	N/A	X336 \$195 ☐	X336 \$195 ☐
	CHOOSE Small trunk support	X327 \$360 ☐	X327 \$360 ☐	N/A
	Large trunk support	N/A	X337 \$415 ☐	X337 \$415 ☐
* Seat	CHOOSE Small	X324 \$330 ☐	N/A	N/A
	Medium	X325 \$385 ☐	X334 \$385 ☐	N/A
	Large	N/A	X335 \$415 ☐	(included)
Headrest (requires trunk support)	CHOOSE Contoured	X376 \$275 ☐	X376 \$275 ☐	X376 \$275 ☐
	Winged	X378 \$300 ☐	X378 \$300 ☐	X378 \$300 ☐
Laterals (pair) (require trunk support)	CHOOSE For small trunk support	X328 \$240 ☐	X328 \$240 ☐	N/A
	For large trunk support	N/A	X338 \$250 ☐	X338 \$250 ☐
Rear handle (rear steer must be factory installed)	CHOOSE Push handle	X312 \$115 ☐	X312 \$115 ☐	X313 \$115 ☐
	Push handle w/ caregiver brake	X314 \$195 ☐	X314 \$195 ☐	X315 \$195 ☐
	Rear steer	X352 \$385 ☐	X353 \$383 ☐	X354 \$383 ☐
	Rear steer w/ caregiver brake	X362 \$470 ☐	X363 \$470 ☐	X364 \$470 ☐
Hand anchors (pair)	CHOOSE Small	X317 \$140 ☐	X317 \$140 ☐	N/A
	Large	N/A	X318 \$140 ☐	X318 \$140 ☐
Front brake (factory installation recommended)	CHOOSE Right hand	X329 \$170 ☐	X339 \$170 ☐	X349 \$170 ☐
	Left hand	X369 \$170 ☐	X379 \$170 ☐	X389 \$170 ☐
Front guide bar		X319 \$195 ☐	X319 \$195 ☐	X319 \$195 ☐
Abductor		X355 \$230 ☐	X356 \$240 ☐	X357 \$250 ☐
Abductor with adduction straps		X365 \$310 ☐	X366 \$320 ☐	X367 \$330 ☐
Pelvic harness	CHOOSE Small	X372 \$200 ☐	N/A	N/A ☐
	Medium	N/A	X373 \$210 ☐	X373 \$210 ☐
	Large	N/A	X374 \$225 ☐	X374 \$225 ☐
Pedal raise kit		X344 \$57 ☐	X344 \$57 ☐	X344 \$57 ☐
Stationary stand		X345 \$115 ☐	X345 \$115 ☐	X345 \$115 ☐
Mounting bar (for device holder)		X346 \$57 ☐	X346 \$57 ☐	X346 \$57 ☐
Front basket		X347 \$57 ☐	X347 \$57 ☐	X347 \$57 ☐
Outdoor cover		X348 \$385 ☐	X358 \$415 ☐	X368 \$415 ☐

* Required



Scan for sizing and pricing:
<https://www.rifton.com/trike>



Small and medium Trikes come in your choice of five colors.




Blue



Red



Teal



Lime



Pink

Large Trikes come in your choice of five colors.



Blue



Red



Teal



White



Black

Tricycle dimensions

User dimensions (inches)	X320 small	X330 medium	X340 large
Inside leg length	12 – 22	19 – 28	25 – 38
 Key user dimension: inside leg length The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place. Important: User's weight must not exceed the maximum working load. Rifton Tricycles are not intended for clients with strong reflexes or poorly developed trunk balance and head control. Consult your client's therapist before ordering.			
Item dimensions (inches)	X320 small	X330 medium	X340 large
Overall width	25½	27½	31½
Overall length	44	53	65
Overall height (at lowest position with wheels / front end removed for transport)	21¼	21¼	30
Overall length with no rear handle (wheels and front end removed)	36	41	49
Overall length with rear handle (wheels and front end removed)	48	51	60
Wheel diameter	12	15	18
Seat center to extended pedal	12 – 22	19 – 28	25 – 36
Seat to top of trunk support (max. height)	16	24	24
Gear ratio	1:1	1:1	1:1
Item weight (lb)	47	56	65
Adjustment range of laterals	small trunk support: 6 – 12 small trunk support: 6 – 12 large trunk support: 9½ – 15½		
Max. working load (lb)	125	250	310