



Nitro HD Rollator

Item #	Color	UOM
10266HD-BR	Brown	1/ea

- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button.
- Back support height easily adjusts with tool-free thumb screw.
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage.
- Plastic clip holds the Nitro closed when frame is folded.
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded.
- Lightweight, aluminum frame.

specifications

Handle Height	36" - 41"
Unit Dimensions	28" (L) x 28.25" (W)
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Weight	28.5 lbs.
Weight Limit	450 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake Handle & Housing)	5 Years