



**Outlander  
Exterior Lift Series**

Order Form / Price List



**Pride Mobility Products Corporation**

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

Account Number: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

Date: \_\_\_\_\_

PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Terms: \_\_\_\_\_

Freight: \_\_\_\_\_

Taken By: \_\_\_\_\_

Vehicle Make: \_\_\_\_\_ Model: \_\_\_\_\_ Year: \_\_\_\_\_

Please read all footnotes - May change vehicle application

**CHOOSE LIFT, OPTIONS & ACCESSORIES**

MODEL NAME	MSRP	WEIGHT(A)	WEIGHT(B)				HITCH TUBE	POWER OPTIONS		
↓	↓	↓	→	8lbs.	9lbs.	11lbs.	9lbs.	N/A	17lbs.	17lbs.
			CLASS II MSRP = \$0	CLASS III/IV MSRP = \$0	CLASS III/IV EXTENDER MSRP = \$213	UNIVERSAL RISER MSRP = \$89	VEHICLE HARD WIRE KIT <sup>6</sup> MSRP = \$0	COMBINATION HARD. PACK w/ BATTERY MSRP = \$330	BATTERY PACK w/ CHARGER <sup>4</sup> MSRP = \$320	
<input type="checkbox"/> OUTLANDER	\$2,050	93lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER DE	\$2,050	88lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER XL	\$2,050	103lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER LP 10"	\$2,050	74lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER LP 14"	\$2,050	73lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER TM <sup>6,7,8</sup>	\$1,595	76lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

SECUREMENT ACCESSORIES - (CHOOSE ONE)	MSRP	WEIGHT(C)	ADDITIONAL SECUREMENT ACCESSORIES	MSRP	WEIGHT(C)
<input type="checkbox"/> Scooter Lockdown Arm (LFTASMB1245)	N/C	8lbs.	<input type="checkbox"/> Scooter Lockdown Arm	\$250	8lbs.
<input type="checkbox"/> 3 Ratcheting Straps	N/C	14lbs.	<input type="checkbox"/> 3 Ratcheting Straps	\$250	14lbs.
<input type="checkbox"/> Powerchair Hold Down Arm (LFTASMB1273) <sup>5</sup>	N/C	6lbs.	<input type="checkbox"/> Powerchair Hold Down Arm (LFTASMB1273) <sup>5</sup>	\$250	6lbs.

**User Information - Vehicle Information**

Hitch to Ground Height: \_\_\_\_\_

**Please read all footnotes - May change vehicle application**

<b>ACCESSORIES -</b>	<b>MSRP</b>	<b>WEIGHT(D)</b>
<input type="checkbox"/> Low Profile Hitch Tube Class III <b>(HTCASMB1022)</b>	\$99.45	9lbs.
<input type="checkbox"/> Hitch Lowering Tube	\$90	11lbs.
<input type="checkbox"/> Swing Away Adaptor <sup>2</sup>	\$600	35lbs.
<input type="checkbox"/> Black, Micro Scooter Weather Cover 61" L x 41.25" H	\$114	5lbs.
<input type="checkbox"/> Black, Medium Scooter Weather Cover 75.75" L x 46.5" H	\$158	6lbs.
<input type="checkbox"/> Black, Large Scooter Weather Cover 84" L x 49" H	\$202	7lbs.
<input type="checkbox"/> Black, Micro Power Chair Weather Cover 62" L x 46.25" H	\$114	5lbs.
<input type="checkbox"/> Black, Medium Power Chair Weather Cover 67" L x 51" H	\$158	6lbs.
<input type="checkbox"/> Black, Large Power Chair Weather Cover 76" L x 52" H	\$202	7lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube <b>(Available as SOF) <sup>1</sup></b>	\$130	5lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 2" Rise <b>(Available as SOF) <sup>1</sup></b>	\$130	6lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 4" Rise <b>(Available as SOF) <sup>1</sup></b>	\$130	6lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 6" Rise <b>(Available as SOF) <sup>1</sup></b>	\$130	7lbs.

**Footnotes:**

1. 10 business day lead time applies.
2. For use when rear door or hatch needs to be opened, swings the exterior lift to one side to allow access w/o removing the lift.
3. Tongue Weight - "weight carrying rating only" - Do not use "weight distribution rating."
4. Must not exceed 380lbs. which is the max weight limit of the lift or the vehicle tongue weight rating (see footnote 3 above).
5. Only option available for the Outlander LP models.
6. Scooter lockdown arm included with the Outlander TM model.
7. Only available with hardwire kit.
8. Low profile Class II hitch included.

**CONFIRM VEHICLE APPLICATION**

° Weights A \_\_\_\_\_ + B \_\_\_\_\_ + C \_\_\_\_\_ + D \_\_\_\_\_ = E \_\_\_\_\_

° E \_\_\_\_\_ + Mobility Device (including all accessory weights) \_\_\_\_\_ = F \_\_\_\_\_

° F \_\_\_\_\_ - Vehicle Tongue Weight<sup>3</sup> \_\_\_\_\_ = G<sup>5</sup> \_\_\_\_\_

**(G must be less than Vehicle Tongue Weight<sup>3</sup> for application to work properly)**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.



OUTLANDER\_LIFT\_ORD\_FORM-07-02-14