

INVACARE®

# Matrix® Personal Back™ 10 Back

## Personal Back™ 10 Plus Back

The Invacare® Matrix® Personal Back™10 and Personal Back™10 Plus Backs combine simplicity of design with effectiveness of application. The highly popular Invacare Personal Back line continues to provide the comfort, postural support and spinal alignment that you have come to expect.

HCPCS Code  
E2615/E2616



Personal Back10 Back



Tall

Regular

Personal Back10 Plus Back

### Invacare® Matrix® Personal Back™10 and Invacare® Matrix® Personal Back™10 Plus Backs

Model no. PBR, PBT

Model no. PPR, PPT

#### Features and Functional Benefits

##### Outer cover:

- Stretchable and breathable cover to allow additional foam inserts

##### Soft foam:

- Foam lining held in place with hook and loop fasteners provides comfort and support
- Allows greater flexibility for addition of foam
- Visco foam upgrade available

##### Back shell:

- Regular and tall heights to accommodate different user postural requirements
- Lightweight aluminum shell that mounts quickly and easily to the wheelchair
- Pre-contoured shell provides excellent contact, improved lateral stability and support
- Deep lateral contour provides increased midline support without interfering with upper extremity range of motion
- Pre-drilled headrest pattern allows mounting of standard headrest mounting hardware

##### Sacral/lumbar support

- Includes sacral support that is easily adjustable to meet individual needs



SACRAL/  
LUMBAR  
SUPPORT

##### Hardware

- Invacare® Quick Set™ Hardware fits canes 7/8 or 1" in diameter
- Can be easily mounted in front of or behind the canes
- Quick release levers for easy removal
- Allows increased adjustability in height, depth (2" in front or behind the canes) and back angle (up to 24°)

#### STANDARD

Width	Regular	Tall
13" – 14"	14"	16"
15" – 26"	16"	19"
Weight capacity		
13" – 20"	250 lb.	
22" – 26"	500 lb.	

#### PLUS

Weight capacity		
13" – 20"	250 lb.	
22" – 26"	500 lb.	



#### FOAM-IN-PLACE KIT

Individually molded insert produced in minutes with user seated in back  
Optional high stretch cover available