



HCPCS
E0570



18080



Available for Private Label
Ask Your Local Drive
Sales Representative

Power Neb Ultra Nebulizer

Item #	Product Description	UOM
18080	Power Neb Ultra Nebulizer with Disposable Neb Kit	1/ea
18081	Power Neb Ultra Nebulizer with Reusable and Disposable Neb Kit	1/ea

- Drive's smallest compressor neb weighs in at a mere 3.2 lbs. and provides a competitive edge price point for providers.
- Built-in neb kit holder for hands-free use.
- Delivers consistent particle size of <5 microns and fast treatment times.
- Complete with disposable neb kit.
- Powerful piston pump is efficient, dependable and quiet.

specifications

Method of Nebulization	Piston Compressor
Particle Size	0.5 µm to 5 µm
Unit Dimensions	6.3" (W) x 3.0" (H) x 5.3" (D)
Weight	3.2 lbs.
Maximum Pressure	30 psi
Operating Pressure	≥ 10 psi
Liter Flow	≥ 8 lpm
Electrical Requirement	115 VAC 60Hz
Warranty	Limited 5 Year

accessories

18090F	Replacement Filter for 18070 and 18080 Character Nebs	10/bag
NEB KIT 500	Disposable Neb Kit with Mouthpiece, "T," and 7' Tubing	50/cs
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs
MASK 001A	Adult Aerosol Mask	50/cs
18031	Universal Nebulizer Carry Shoulder Bag	1/ea