

Jazzy® 1103 ULTRA

The Jazzy® 1103 Ultra brings sporty style and high performance features, like a standard power elevating seat and Active-Trac® Suspension, together into one perfectly balanced power chair. Yet, even fully loaded, the Jazzy 1103 Ultra remains highly compact and maneuverable.



Jazzy® 1103 Ultra shown
in Candy Apple Red with
standard power elevating
gray contoured medium-
back seat, foot platform and
VSI controller.



INNOVATION

Standard power elevating seat provides greater access to your environment



PERFORMANCE

Active-Trac® Suspension for enhanced outdoor performance



STYLE

Sleek styling for a sporty edge



MANEUVERABILITY

Compact design and tight turning radius perfect for small spaces

SPECIFICATIONS

Weight capacity:	300 lbs.		
Maximum speed:	Up to 4 mph		
Ground clearance:	2.5"		
Turning radius:	18.5"		
Overall length:	35.5"		
Overall width:	23"		
Wheels:	Drive: 10" pneumatic Casters: 8" solid Anti-tips: 5" solid		
Suspension:	Active-Trac®		
Drive train:	(2) motor, mid-wheel drive		
Braking system:	Intelligent braking (electronic, regenerative disc brakes)		
Std. Electronics:	50 amp, 6-key VSI		
Specialty controls:	Yes		
Battery charger:	On-board, 4 amp		
Per-charge range:	Up to 25 miles		
Battery requirements:	(2) 12 volt deep cycle		
Battery size:	U-1		
Transport weight:			
Batteries:	24.5 lbs. each		
Base:	107 lbs.		
Standard seat:	33 lbs. (medium-back)		
Std. Seat dimensions:	18" x 18"		
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors		
Colors:			
Candy Apple Red	Viper Blue	Onyx Black	Champagne

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- 300 lbs. weight capacity
- Standard power elevating seat provides greater access to objects out of reach
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- Standard lap belt
- On-board battery charger

OPTIONS

- Contoured high-back seat
- Synergy® Seating System
- Solid seat pan
- Solid tires
- TRU-Balance® Power Positioning System*
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- Saddlebag
- Off-board charger

* Requires off-board charger;
weight capacity is 200 lbs.