

Model	SC7132
Weight Capacity	400 lbs.
Maximum Speed <sup>2</sup>	Up to 9.3 mph
Ground Clearance <sup>4</sup>	3.8" at transaxle
Turning Radius <sup>4</sup>	69.3"
Overall Length <sup>4</sup>	57.5"
Overall Width <sup>4</sup>	29.7"
Seat-to-Ground Height	21.9" - 24.9"
Seat-to-Deck Height	17.9" - 20.9"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Front and rear (CTS Suspension)
Range Per Charge <sup>2,3,7</sup>	Up to 40 miles with upgradable dual lithium battery configuration
Total Weight <sup>6</sup>	With standard battery: 242.6 lbs. Without batteries: 220.6 lbs.
Heaviest Piece When Disassembled <sup>4</sup>	N/A (Unit does not disassemble)
Standard Seating <sup>4</sup>	Type: Limited reclining high back with sliders Weight: 54.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>5</sup>	25 volt, 50Ah lithium ion
Battery Charger	8 amp
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 3 years

### HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## INSERT YOURSELF: ANYWHERE

FDA Class II  
Medical Device<sup>1</sup>



FDA Class II  
Medical Device<sup>1</sup>



**PRIDE**<sup>®</sup>  
M O B I L I T Y

800-800-8586 | pridemobility.com



1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) +10%, -5%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Includes standard seat weight.
6. Battery weight may vary based on battery manufacturer.
7. With dual battery configuration.

# PURSUIT<sup>2</sup>

Engineered for Adventure



40 MILE  
RANGE



LONGER  
LIFE



CRUISING  
SPEED



WEIGHT  
CAPACITY

LIVE YOUR BEST<sup>®</sup>

pridemobility.com |    

FDA Class II Medical Device<sup>1</sup>

# PURSUIT<sup>2</sup>

Built to accompany you on all your adventures, big or small.

Memory foam seat with 6" slider for adjustable support and premium, all-day comfort

Durable plastic shrouds

Tiller adjustment lever and USB charger

Front basket to hold whatever you need for the road ahead

Enhanced LED lighting

Pneumatic tires mounted on split rims for off-road thrills


FDA Class II Medical Device<sup>1</sup>



Infinitely adjustable tiller with USB charger



Enhanced LED lighting



Front and rear CTS suspension

## ENGINEERED FOR ADVENTURE

The Pursuit 2 proves the best adventures are yet to come as its fast cruising speed, weight capacity and pneumatic tires can tackle any road (or backroad).

## FEATURES

- 400 lb. weight capacity
- 9.3 mph top cruising speed
- Up to a 40-mile range with upgradable dual lithium battery configuration
- Durable plastic shrouds
- Memory foam seat with 6" slider for adjustable support and premium, all-day comfort
- Infinitely adjustable tiller with USB charger
- Enhanced LED lighting for clear visibility
- Pneumatic tires for varied terrains and surfaces
- 3-year warranty on lithium battery
- Depth adjustable seat options up to 22" x 22"

## AVAILABLE COLORS

Black

Gray (Matte)

Orange

## ACCESSORIES\*\*

- Wishbone Crutch Holder
- 3" Single/Double Cane or Crutch Holder
- Cane Crutch Holder
- Cell Phone Holder
- Cup Holder\*
- Forearm Crutch Holder
- Lap Belt
- Military Patches
- Oxygen Tank Holder
- Rear Basket\*
- Saddlebag\*
- Safety Flag
- Storage Pod
- Walker Holder
- Weather Cover

**NEW**

Storage pod attaches to the rear of seat.



\* Cup holder, rear basket and saddlebag can be ordered as a bundle.

\*\* Only one rear accessory can be used at a time.