

The Power Within



portable power **QUICKIE® MELODY™**

QUICKIE®
POWER CHAIRS

Fast and very compact, the Quickie® Melody™ quickly disassembles into four easy to manage components, each weighing less than 36 pounds! As a lightweight, highly portable power chair, the Melody sets a new standard for versatility that on-the-go people will love.

Quickie Melody Power Wheelchair

Extremely light, modular and manageable

Modular Components



It's a remarkable revolution in portability and power: the Quickie Melody. This lightweight chair **quickly disassembles** into four components (two base sections, seat and battery) for **easy transport**. And no piece weighs over 36 pounds. Yet it's got the power and speed you need, up to four mph, plus a turning radius of just 24.5 inches.

It's **stylish and compact**, with a small footprint.


Quick-Release Battery



For even greater **ease when traveling**, you can charge the Melody through the joystick or battery port. There's also a **quick-release, wireless battery** for easier access that eliminates the hassle of undoing battery cables prior to transport.

For an ideal fit, there's a two-inch center-of-gravity adjustment, two seat-to-floor heights (21 and 22 inches), and an **easy-release, swivel seat** that locks in four directions for safe access. The Quickie Melody is a perfect choice for people who love to travel and Get Out There!

To learn more, contact your local dealer, call Sunrise Medical at **1-800-333-4000** or visit **www.sunrisemedical.com**.



portable power **QUICKIE® MELODY™**

QUICKIE®
POWER CHAIRS

Melody

Speed	4 mph
Turning Radius	24.5"
Ground Clearance	2"
Maximum Range	Up to 10 Miles
Chair Length (w/footplate)	32"
Chair Width	22.5"
Seat-to-Floor Height	21" and 22"
Overall Weight	113 lbs (90 lbs w/o Battery)
Heaviest Component	35.5 lbs
Weight Capacity	250 lbs
Upholstery	Black Vinyl
Seat Type	Low Back Comfort Seat
Drive Tire Size	8" Solid
Front Caster/Wheel Size	6" Solid
Shroud Color	Red
Battery Type	2 ea 12 Ahr
Charger Location	Off-board
Charger Output	2 Amp
Controller Type	PG VR2
Controller Output	60 Amp

WARRANTY

Base Frame	5 Years
Controller	2 Years
Motors, Gear Box and Remote Joystick	1 Year
Wheels, Upholstery, Plastic or Rubber Parts and Painted Surfaces	3 Months

