

Premium Series Shower Chair with Back and Arms

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easy, tool-free height adjustment with height indicator
- Adjusts in ½" increments
- Enhanced comfort seat fits all shapes and sizes
- Easy, tool-free removable back, legs and arms
- Lightweight design has large stable legs with soft non-slip feet
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
RTL12505	BATH SEAT WITH BACK, 1ea, 4/cv
SPECIFICATIONS	Description
Warranty	Limited lifetime
Seat Dimensions	21.75" (W) x 16.5" (D)
Weight Capacity	350 lbs
Carton Shipping Weight	48 lbs.
Weight	11.5 lb