

SILVER SPORT FULL RECLINING



**SILVER VEIN
FINISH ON
FRAME.**

Silver Sport Full Reclining Wheelchair Single Axle (K0007)

FEATURES:

- New state-of-the-art hydraulic reclining mechanism allows for infinite adjustments up to 180°.
- Arms have support extension for reclining position.
- Carbon Steel frame with silver vein finish for an attractive, maintainable finish.
- Durable, lightweight nylon upholstery.
- Padded armrests provide added patient comfort.
- Composite Mag-style wheels are lightweight and maintenance free.
- Comes standard with swing-away elevating legrests.
- Wheels set back on frame prevents tipping.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Plastic footplates are attractive, crack-proof and lightweight.
- Rear anti-tippers standard.
- 8" front casters.
- Headrest extension with cushioned head immobilizer standard.
- Comes with push to lock wheel locks.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.