2024 order form for small

Use dimension chart to select appropriate size.



Elf-6-789 Call: 800.571.8198			Standard R820 small	A STATE OF THE PARTY OF THE PAR	J	Hi/lo R830 small		d
*Seat and back includes seatbelt	CHO	Adjustable	R821	\$650		R831	\$1,330	
(Hi/lo includes footboard and push handles)	CHOOSE	Adjustable w/ back spring	R822	\$745		R832	\$1,450	
*Base	E	Without spring	R823	\$380		R833	\$2,270	
	CHOOSE	With spring	R824	\$450		N/A		
*Arm supports (pair)	운	Armrests	R827	\$265		R827	\$265	
	CHOOSE	Forearm prompts	R809	\$650		R809	\$650	
*Adjustable legs		Short	R803	\$55		N/A		
	윤	Long	R894	\$85		N/A	•••••	
	CHOOSE	Short with casters	R805	\$325		N/A	•••••	
		Long with casters	R896	\$360		N/A		
* Pads		Tan	R815	\$270		R815	\$270	
		Blue	R816	\$270		R816	\$270	
	CHOOSE	Red	R817	\$270		R817	\$270	
	OSE	Pink	R971	\$270		R971	\$270	
		Purple	R972	\$270		R972	\$270	
		Green	R973	\$270		R973	\$270	
Headrest	0	Contoured	R897	\$260		R897	\$260	
	CHOOSE	Flat	R898	\$260		R898	\$260	
	m	Adjustable winged	R881	\$290		R881	\$290	
Laterals (pair)	CHOOSE	Small	R828	\$220		R828	\$220	
	3SOC	Small with chest strap	R819	\$265		R819	\$265	
Push handles (recommended for all chairs with wheels)			R878	\$166		Included		
Tray (requires armrests)			R667	\$265		R667	\$265	
Handhold			R649	\$31		R649	\$31	
Chest strap (wide)			R802	\$95		R802	\$95	
Butterfly harness			R834	\$130		R834	\$130	
Thigh belt			R811	\$100		R811	\$100	
Pelvic harness			R845	\$170		R845	\$170	
Hip guides (pair)			R814	\$225		R814	\$225	
Abductor			R884	\$140		R884	\$140	
Adductors (pair)			R892	\$220		R892	\$220	
Leg prompt			R854	\$235		R854	\$235	
Footboard (recommended for chairs with long legs)			R847	\$515		Included		
Mini kit (backrest insert and trunk support)			R812	\$140		R812	\$140	
Ankle straps (pair)			R899	\$100		R899	\$100	
Sandals (pair)			R944	\$100		R944	\$100	
Wedges (pair, require sandals)			R697	\$74		R697	\$74	
Lumbar and seat support kit			R886	\$73		R886	\$73	
Backrest filler pad			R849	\$135		R849	\$135	
Whitmyer headrest adapter			R804	\$63		R804	\$63	

^{*} Required



Room to grow

Small chair with mini kit at its smallest dimensions

- Backrest height 12½"
- Trunk support width 5½"
- Seat depth 7"
- 4 Seat width 7"
- Armrest height above seat 5"
- Seat height above footboard 7"



Small chair at its largest dimensions

- Backrest height 15½"
- 2 Distance between laterals 11½"
- Seat depth 12"
- Seat width 12"
- ⑤ Armrest height above seat 7½"
- Seat height above footboard 12"





Small standard base

R820 Rifton Activity Chair

User dimensions (inches)	
Height with mini kit	32-48 28-38



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)	R820 small Standard base
Frame width	short legs: 21
	long legs: 23
	short legs w/ casters: 221/2
	long legs w/ casters: 231/2
Seat height above floor	short legs: 9½-12½
	long legs: 18½-21½
	short legs w/ casters: 13½-16½
	long legs w/ casters: 18½-21½
Seat angle (tilt-in-space) – degrees	15° forward, 15° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	7-12
Seat width with hip guides without hip guides	7–9 12
Seat depth with mini kit	8-12 7-11
Armrest height above seat	5-71/2
Trunk support width	5½-11½
Backrest height	121/2 – 151/2
Headrest height above seat	14½-21
Max. working load (lb)	75



Small hi/lo base

R830 Rifton Activity Chair

U	ser dimensions (inches)	
He	eight ith mini kit	32-48 28-38



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)	R830 small Hi/lo base
Frame width	26
Seat height above floor	10-25
Seat angle (tilt-in-space) – degrees	15° forward, 25° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45°-110°
Seat height above footboard	7-12
Seat width with hip guides without hip guides	7-9 12
Seat depth with mini kit	8-12 7-11
Armrest height above seat	5-7½
Trunk support width	5½-11½
Backrest height	12½-15½
Headrest height above seat	141/2-21
Max. working load (lb)	75

