

2026 order form

Use dimension chart to select appropriate size.

Scan for sizing and pricing:
<https://www.rifton.com/standers>



		S410 size 1	S420 size 2	S430 size 3
* Frame		S411 \$1,990 ☐	S421 \$2,640 ☐	S431 \$3,270 ☐
* Color		Blue ☐	Blue ☐	Blue ☐
		Red ☐	Red ☐	Red ☐
		Lime ☐	N/A	N/A
* CHOOSE CONFIGURATION	Prone (with pelvic harness)	S412 \$315 ☐	S422 \$325 ☐	S432 \$335 ☐
	New! Prone upper trunk support (recommended)	N/A	N/A	S408 \$260 ☐
	* Trunk strap ☐	Standard	S462 \$315 ☐	S463 \$325 ☐
		Rotation control	S466 \$315 ☐	S467 \$325 ☐
	Prone knee cuffs (pair)	S485 \$210 ☐	S485 \$210 ☐	S486 \$235 ☐
	Supine (with upper trunk support)	S413 \$370 ☐	S423 \$380 ☐	S433 \$390 ☐
	Head laterals (pair) ☐	Small	S497 \$210 ☐	S497 \$210 ☐
		Large	N/A	S498 \$225 ☐
	* Trunk straps ☐	Butterfly harness (standard)	S456 \$320 ☐	S457 \$330 ☐
		Butterfly harness (slim-cut)	N/A	S459 \$340 ☐
		Trunk strap (standard)	S462 \$315 ☐	S463 \$325 ☐
		Trunk strap (rotation control)	S466 \$315 ☐	S467 \$325 ☐
	* Hip strap ☐	Standard	S472 \$315 ☐	S473 \$325 ☐
		Rotation control	S476 \$315 ☐	S477 \$325 ☐
	* Knee supports ☐	Soft straps (pair)	S401 \$210 ☐	S402 \$225 ☐
		Split supports (pair)	N/A	S487 \$450 ☐
	Multi-Position (with head support, pelvic harness)	S414 \$690 ☐	S424 \$710 ☐	S434 \$730 ☐
	Head laterals ☐	Small	S497 \$210 ☐	S497 \$210 ☐
		Large	N/A	S498 \$225 ☐
	Butterfly harness ☐	Standard	S456 \$320 ☐	S457 \$330 ☐
		Slim-cut	N/A	S459 \$340 ☐
	* Trunk strap ☐	Standard	S462 \$315 ☐	S463 \$325 ☐
		Rotation control	S466 \$315 ☐	S467 \$325 ☐
	* Hip strap ☐	Standard	S472 \$315 ☐	S473 \$325 ☐
	Rotation control	S476 \$315 ☐	S477 \$325 ☐	
* Knee supports ☐	Soft straps (pair)	S401 \$210 ☐	S402 \$225 ☐	
	Split supports (pair)	N/A	S487 \$450 ☐	
	New! Prone upper trunk support (recommended)	N/A	N/A	S408 \$260 ☐
	Prone knee cuffs (pair)	S485 \$210 ☐	S485 \$210 ☐	S486 \$235 ☐
* Trunk laterals ☐	Small	S494 \$210 ☐	S494 \$210 ☐	N/A
	Large	N/A	S495 \$225 ☐	S495 \$225 ☐
* Leg type ☐	Standard	S415 \$700 ☐	S425 \$720 ☐	S435 \$760 ☐
	Advanced	N/A	S426 \$1,030 ☐	S436 \$1,090 ☐
Tray		S491 \$580 ☐	S491 \$580 ☐	S492 \$630 ☐
Additional trunk laterals ☐	Small	S494 \$210 ☐	S494 \$210 ☐	N/A
	Large	N/A	S495 \$225 ☐	S495 \$225 ☐
Medial/lateral thigh supports (pair) ☐	Small	S482 \$210 ☐	S482 \$210 ☐	N/A
	Large	N/A	S483 \$225 ☐	S483 \$225 ☐
Additional medial/lateral thigh supports (pair) ☐	Small	S482 \$210 ☐	S482 \$210 ☐	N/A
	Large	N/A	S483 \$225 ☐	S483 \$225 ☐
Size 1 sandal lift kit		S409 \$57 ☐	NA	NA
Headrest adapter plate kit		N/A	S427 \$265 ☐	S427 \$265 ☐

* Required

Stander dimensions

User dimensions (inches)	S410 Size 1	S420 Size 2	S430 Size 3
Height	28–42*	38–53	48–65
 Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth. *The smallest users may require the S409 Sandal lift kit for best positioning. Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)	S410 Size 1	S420 Size 2	S430 Size 3
Base length x width	32½ x 21½	39 x 27	43 x 29
Height of horizontal body support above floor	22–28	22½–36½	23½–38
Sandal to knee height	6¾–11	10–15¼	12¼–19
Knee to hip height	6¼–10½	8½–13½	12–17¼
Sandal to top of headrest	31½–43	41–55	49–66½
Width between trunk laterals	5½–10	7–11½	9–14
Width between medial/lateral thigh supports	2½–4	3¼–5	4–5½
Basic item weight PRONE (lb)	47	72	104
Basic item weight SUPINE (lb)	52	78	115
Max. working load (lb)	75	100	185



New! Size 3 Stander at its smallest and largest settings